

# **MICROLAB® 1000**

## **INSTRUCTION MANUAL**

2988

## TABLE OF CONTENTS

|                                                   | <u>PAGE</u> |
|---------------------------------------------------|-------------|
| Forword                                           | 1           |
| Warranty                                          | 2           |
| Warranty Registration Form                        | 3           |
| Description and location of components            | 4           |
| Unpacking/Electrical installation                 | 7           |
| Introduction                                      | 8           |
| Operating instructions (new user)                 | 10          |
| OPERATING PROCEDURE                               |             |
| Symbols used                                      | 14          |
| Program as diluter                                | 15          |
| Program as double diluter                         | 20          |
| Program as single or multiple dispenser           | 25          |
| Program as a serum distributor                    | 29          |
| Program as serial diluter                         | 35          |
| Program as serial dispenser                       | 41          |
| Program as multiple pick up                       | 46          |
| To find a stored method or an unoccupied location | 51          |
| To delete one or all programmes                   | 52          |
|                                                   |             |
| Fitting the valve block                           | 53          |
| Fitting the syringes                              | 53          |
| To exchange the valve block                       | 53          |
| To exchange syringes                              | 54          |
| To replace defective syringes                     | 54          |
| Aspiration tube                                   | 54          |
| Sample/delivery tube                              | 54          |
| Selection of tubing                               | 55          |
| Cleaning of the FEP tubing                        | 55          |
| Technical specifications                          | 56          |
| Figure 9.1.1 (front view)                         | 58          |
| Figure 9.1.2 (rear view)                          | 59          |
| Description of the RS 232 C interface             | 60          |
| Test certificate                                  | 67          |

## SECTION 1

### Forword

- 1.1 Precision and accuracy in measuring microlitre volumes are difficult to obtain in a convenient manner. That is why HAMILTON produced the MICROLAB 1000. Its performance and simplicity of operation when used for diluting, pipetting, or distributing samples and reagents brings microlitre liquid handling problems down to earth.

With the MICROLAB 1000 complete methods may be stored to be recalled at anytime. The built in microprocessor will carry out the task with the same high degree of precision and accuracy expected from HAMILTON products, bringing a new dimension to sample processing.

HAMILTON have been manufacturing microlitre syringes for nearly 30 years and have built up a world wide reputation with literally millions of products in use including much instrumentation. The MICROLAB 1000 is one of a series of instruments culminating from this vast experience. May we congratulate you on becoming the owner of this most up to date, microprocessor controlled liquid handling system.

## **Warranty**

- 1.2 HAMILTON warrants all instruments to be free from any defects in construction, material or workmanship for a period of 18 months from the date of shipment ex factory.

HAMILTON will repair any defective instrument or part thereof if the device is returned within the warranty period to the factory or an authorized dealer. Seller shall not be liable nor responsible for any incidental or consequential damages.

HAMILTON shall be relieved of any liability under this warranty if the instrument is not used in accordance with manufacturer's instructions, not regularly maintained or used for a purpose for which it was not designed. Any product which has been mechanically or electronically altered without specific authorization is excluded from this warranty.

Expendable items such as syringes, tubing and valves are warranted at the time of delivery only.

Please return defective instruments carefully packed in the original boxes to the authorized dealer or to the factory. Freight charges are at buyer's expense.

For any assistance contact your nearest HAMILTON dealer.

Manufactured by Hamilton Bonaduz AG, CH-7402 Bonaduz, Switzerland

ML 1000/10.84.150  
610 047



## Warranty Registration Form

Please fill in this form carefully and mail it to the address listed below:

MICROLAB 1000 Serial No. \_\_\_\_\_

Date of installation: \_\_\_\_\_

Purchased from: \_\_\_\_\_

My specific application is following: \_\_\_\_\_

Name: \_\_\_\_\_

Address: \_\_\_\_\_

Country: \_\_\_\_\_

Mailing address:

Hamilton Bonaduz AG  
P.O. Box 26

CH-7402 Bonaduz  
Switzerland

## SECTION 2

### Description and location of components

Refer to figures 9.1.1 and 9.1.2

- 2.1 Fuses  
The instrument is equipped with two fuses, located in the main socket, one in the line and the other in the neutral. In case of a short circuit, the fuse in the line will blow first. The line fuse is located in the mains filter unit and is accessible by pulling down mains adjustment panel. With mains connector removed insert a small screwdriver in slot and pull panel out horizontally. The fuse holder assembly is indicated with an arrow. Remove the assembly by pulling out the unit. After inserting the fuse make certain that the arrow marked on the holder assembly is directed according to the indication on the cover. Fuse rating for 220-240 volt operation is 1 amp and for 100-120 volt 2 amp.
- 2.2 Mains voltage adjust.  
Situating on the mains filter unit at the back panel, four nominal voltages can be selected. 100-120-220-240 input voltage must not exceed  $\pm 10\%$  of nominal. To change operating voltage pull out the panel using a small screwdriver in the slot above mains socket and remove the voltage selector drum, turn it to the desired mains voltage. After closing the panel, the selected voltage is indicated in the window. Fuse ratings change with voltage see under fuses 2.1.
- 2.3 Mains connector  
A 3-pin moulded plug is supplied with each unit, polarised to ensure correct connection to line (phase) and neutral inputs.
- 2.4 Mains switch  
The mains switch is a double pole type controlling both line and neutral inputs. Press the top to switch on. When read circle is visible switch is in the off position.
- 2.5 Display  
20 character, alpha-numeric, fluorescent display, used to perform dialogue between user and microprocessor.
- 2.6 Display-cover  
This carries a special filter for the display. It can be lifted to give access to the valve block and volume control for audible signal.
- 2.7 Touch-board  
Touch type key-board is used to enter numbers and answer dialogue questions.
- 2.8 Valveblock  
Situating under the display-cover. Two types are used depending on volumes used and designated MICRO and MACRO. The built in microprocessor will inform the user which valve should be used, for any application.

**Important:** The valveblock must only be changed when the microprocessor requests it. The type fitted is retained in memory and will be taken into account before each run.

## 2.9 Syringe drives

Are controlled only by the computer and must never be forced manually. The left hand side is designated as A-syringe for reagent/diluent measurement and the right hand side as B-syringe for sample measurement. Each drive has a resolution of syringe total capacity  
1000

**Important:** The syringes must only be changed when the microprocessor requests it. The types fitted are retained in memory and will be taken into account before each run.

## 2.10 Handactuator

Unit containing the delivery tip, ready light and START button. It is used to control directly the aspiration and dispensing of solutions in normal operation. A small lamp indicates when instrument is ready to proceed. The start button is also used during setting-up, to install the syringes, for zero-setting and to control the flush cycle. All operations are indicated on the display. The handactuator is plugged-in in the rear panel socket marked AUX.

## 2.11 Serial interface

RS 232-C

## 2.12 Audio Signal

Acceptance of key-board commands are indicated by sounding a short bleep. Errors are announced by a long bleep with an indication of the cause briefly shown on the display. The volume of this signal may be adjusted using the control located to the right of the display. Clockwise rotation will increase volume.

## 2.13 Syringes

Various combinations of syringes are possible and may be selected by the user, but the microprocessor will, for any application, inform the user of the selection for optimum performance.

**Important:** The syringes must only be changed when the microprocessor requests it. The types fitted are retained in memory and will be taken into account before each run.

## 2.14 Tubing

Standard tubing is supplied, a reagent aspirate tubing 2 x 1 mm dia, 500 mm long and a delivery tubing 2 x 1 mm, 900 mm long and are fitted to the left side of the valveblock and to the sample syringe respectively. If larger samples than 700 ul are contemplated then a 2 x 3 mm, 900 mm long sample tubing should be used. In any event it is not recommended that samples be taken into the sample syringe as unacceptable carry-over may result.

Note

- Capital letters in display indicate operator activity - lower case indicate instrument activity.
- ? on the display indicates yes or no response is required by the user.
- Flashing word indicates option can be selected with yes response or rejected by no. Response, also when start button should be used.
- Green light on handactuator indicates start button is active, extinguishes when command accepted.
- Short audible bleep indicates key-board command accepted.
- Long audible bleep indicates error detected or overload has occurred.

## **SECTION 3**

### **3.1 Unpacking**

The MICROLAB 1000 is a precision instrument, therefore, the packing material has been selected carefully to provide maximum protection during transportation. Upon receipt, the packing material should be examined for indications of damage. In the case of damage occurring during transportation, contact the carrier immediately otherwise contact your HAMILTON dealer.

### **3.2 Electrical installation**

Unpack the MICROLAB 1000 and all accessories from the box and inspect for physical damage. Report any shipping damage immediately.

Place the MICROLAB 1000 in the position where it is to be operated, but not directly over heat sources or in direct sunlight (temperature extremes).

Ensure mains voltage is set correctly, refer to 2.2.

Ensure correct fuse is fitted, refer to 2.1.

Fit handactuator plug into AUX socket on rear panel, see 2.10.

Connect mains supply, cable supplied with polarised input, see 2.3.

Serial interface connection, see 10.1.

## SECTION 4

### 4.1 Introduction

The MICROLAB 1000 has eight basic functions. It can operate as a diluter, double diluter, serial diluter, dispenser, multiple dispenser, serum distributor or in multiple pick up mode. In the diluter mode it will deliver diluted samples in a pre-determined ratio using one diluent/reagent, as a double diluter by using two diluents/reagents, and as serial diluter selecting up to ten different dilution ratios. In the dispenser mode the instrument will accurately pipette aliquots of reagent from a stock solution. When used as a serum distributor it will deliver the sample into up to ten receptacles, each of different volume if required, and in multiple pick up mode up to ten different samples can be aspirated and dispensed together.

Premature mixing of reagents and sample is avoided by physical separation of liquids by using a programmed air gap, with cross contamination minimised in the serum distributor mode by including a wash after each cycle.

To obtain the highest performance the unit uses two measuring syringes operating in a completely hydraulic system. Two syringes allow large ratio dilutions to be made, up to 1:5000, without loss of performance. The diluent syringe, left hand side, can be large, up to 25000  $\mu$ l and the sample syringe, right hand side, small, thus preserving good resolution for small sample. It is not intended that the sample should ever be drawn into the syringe, but that it be contained within the delivery tubing. The standard tubing, 900 mm in length, will hold up to 700  $\mu$ l. This method of operation practically eliminates carry-over by the washing action of the diluent each time a dilution is made.

All quantities and instructions are entered through the touch-key-board and memorised for future use. This information will be filed under a method number allocated by the user to enable immediate recall at any time. The data will be retained even if the power is removed from the instrument.

Up to fifty methods can be permanently stored at one time, but it is a simple process to delete or amend methods, should the need arise. An overview of each of the stored methods can easily be displayed.

The operator controls the MICROLAB 1000 by answering YES or NO to options offered, which are indicated by causing the option to flash in the display, and by answering questions indicated by a question mark. If the option offered is the one required, it can be selected by using the YES key. The NO key will reject this option and present the next. When the start button on the handactuator should be used, it is indicated in the same way, using a flashing word. When a command from the user is accepted a short bleep will sound but if an error or overload occurs a long bleep will sound together with a brief indication of the cause on the display.

To become familiar with this technique, proceed with the following exercise:

Switch on the MICROLAB 1000 at the rear panel. On the display "SELFTEST COMPLETED" will appear for a few seconds. The instrument is now ready for use. If the instrument has been switched off during programming, the display will indicate "USER MEMORY FAILURE". Press the CE-key to reset the instrument. It should show now **RUN** OR PROG OR AID. The box around a word will be used throughout this text to indicate that this is the flashing word. Touch the NO key, when the command is accepted a short bleep will sound and the display will change to RUN OR **PROG** OR AID. This means NO the RUN channel is not required, so now the microprocessor is asking do you want PROG mode? Touch NO again to change option to RUN OR PROG OR **AID**. Assume this mode is required, so touch YES. Display will change to REVIEW OR DELETE. Touch YES meaning the REVIEW feature is required. Since no methods as yet have been entered the microprocessor will inform the user - location 1 "NO METHOD". Each time YES is touched the location will advance, allowing an overview of the content of each of the fifty locations, so touch YES fifty more times to return to **REVIEW** OR DELETE. Touch NO twice to instruct the microprocessor that they are not wanted and return the dialogue to the starting point: **RUN** OR PROG OR AID.



## SECTION 5

### 5.1 Operating instructions (new user)

In general the MICROLAB 1000 will lead the user through a simple procedure as discussed in the introduction above. This procedure follows a pattern and is easily learned. The operation is divided in three channels. The RUN channel will be used each time the instrument is asked to routinely perform one of the stored programmes. PROG channel will be used to store new methods or to amend the existing ones, and the AID channel allows the user to look at an overview on the display of each method stored, or to delete an entire method, freeing that location for storing a new method.

With many alternative sub-routines possible the easiest way to understand the MICROLAB 1000 is by example. The following example is designed to teach a new user how to program and run the instrument. After working through this exercise it should be possible, using the procedures under section 6, to use the instrument in its other functions. The operator can be confident in operating this unit since any errors will be detected and announced.

When first switched on, the microprocessor will perform a selfcheck, if all is in order the display will briefly inform user - SELF TEST COMPLETED - if an error is found the message will read - USER MEMORY FAILURE. Simply touch any key and the offending method will be deleted from memory.

### 5.2 Programming as a multiple dispenser

#### DISPLAY

- |       |                                                                                                                                                                                  |                                                                                                                                                                                                                                                                                         |
|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 5.2.1 | Switch on mains supply and select PROG mode.<br>Touch NO and then YES (brackets show user activity)                                                                              | <div style="border: 1px solid black; padding: 2px; display: inline-block;">RUN</div> OR <div style="border: 1px solid black; padding: 2px; display: inline-block;">PROG</div> OR AID<br>(NO) (YES)                                                                                      |
| 5.2.2 | Enter a number between 1 and 50 to identify this method for future recall. It should be an unoccupied location, say 20, so touch "2", "0", "YES" keys.                           | METHOD #<br>(2,0,YES)                                                                                                                                                                                                                                                                   |
| 5.2.3 | Select dispenser (DISP) from the four functions offered. DiLuter, Double DiLuter, DISPenser and DISTributor. Touch NO twice and YES once.                                        | <div style="border: 1px solid black; padding: 2px; display: inline-block;">DIL</div> <div style="border: 1px solid black; padding: 2px; display: inline-block;">DDL</div> <div style="border: 1px solid black; padding: 2px; display: inline-block;">DISP</div> DIST<br>(NO) (NO) (YES) |
| 5.2.4 | Select REP.DISP. Select aspiration automatically. The filling of the syringe is executed automatically, or manually, by touching NO. Touch YES.                                  | <div style="border: 1px solid black; padding: 2px; display: inline-block;">REP.DISP</div> / SERIAL DISP.<br>(YES)                                                                                                                                                                       |
| 5.2.5 | Select dispensing automatically. In automatic operation the aliquots will be dispensed after a set time between 0.1 sec. to 9.9 sec. without user indicating each one. Touch YES | (YES)                                                                                                                                                                                                                                                                                   |
| 5.2.6 | Select the time delay between each dispensing. One second is shown as a default. Enter "0" "." "8". A delay of 0.8 sec. is now programmed.                                       | (0, ., 8, YES)                                                                                                                                                                                                                                                                          |



- 5.2.7 Enter the desired size of each aliquot - say 100  $\mu$ L. REAGENT 100  $\mu$ L  
Touch "1", "0", "0", "YES". (1,0,0 YES)
- 5.2.8 Enter number of aliquots required up to 100 - say 10. #OF ALIQUOTS  
Touch "1", "0", "YES". (1,0, YES)
- 5.2.9 MICROLAB 1000 is now programmed to aspirate 1000  $\mu$ L of reagent and dispense it in 10 steps of 100  $\mu$ L each.
- 5.2.10 The microprocessor has calculated the size of the liquid system required for optimum performance and is now giving the user the opportunity to accept or change these values. Select parameter verification by YES. VERIFY?  
(YES)
- 5.2.11 Select syringe verification by YES. Values can be accepted without verification by YES. :FUNCT:VOL:SYR:SPEED  
(NO) (NO)(YES)
- 5.2.12 Display shows optimal computed syringe size of A-SYRINGE AD  
1.0 mL  
1.0 mL, by responding with NO next suitable syringe size will be substituted but some performance will be lost. (NO)
- 5.2.13 Accept this size by YES. The NO key can be pressed again and again until a suitable size for user is shown, use YES to accept. A-SYRINGE AD 2.5 mL  
(YES)
- 5.2.14 Accept with YES. NO will substitute alternative but some performance will be lost. MICRO-VALVE  
(YES)
- 5.2.15 In the same way computed times for syringe movement can be verified or changed. If bubbles form in liquid system during aspiration, the speed should be reduced. Touch YES. FUNCT VOL SYR SPEED  
(YES)
- 5.2.16 The time indicated on the display is the time taken to aspirate the total syringe volume. In this case only 1 mL will be aspirated into a 2.5 mL syringe, so the time taken will be 4 sec.  $\times \frac{1.0 \text{ mL}}{2.5 \text{ mL}} = 1.6 \text{ secs}$ . To increase the speed a shorter time should be selected. Change this parameter by entering "3". Aspirate now will be 3 sec.  $\times \frac{1.0 \text{ mL}}{2.5 \text{ mL}} = 1.2 \text{ secs}$ . asp REAGENT 4 sec  
(NO, 3, YES)
- 5.2.17 Accept dispense speed of 2 sec.  $\times \frac{1.0 \text{ mL}}{2.5 \text{ mL}} = 0.8 \text{ secs}$  by disp REAGENT 2 sec  
(YES)  
touching YES. This could be changed as in 5.2.16 above.
- 5.2.18 The microprocessor is offering the opportunity to check or change the parameters set for this method. VERIFY?  
Touch NO and exit the PROG mode. YES would allow verification or ammendment to be made. (NO)
- 5.2.19 Programming as a multiple dispenser is now completed RUN OR PROG OR AID  
and filed under method number 20.

- 5.3 Running a stored method (in conjunction with 5.2)** DISPLAY
- 5.3.1 Select run channel **RUN** OR PROG OR AID (YES)
- 5.3.2 YES to run method shown and NO to run another method METHOD NO 20  
20 RDISP aut 10x100 (YES)
- 5.3.3 This means syringes already fitted are of the incorrect size and must be replaced. Press START button on handactuator to lower drives to allow syringes to be fitted. If syringes are already fitted and are suitable, programme will continue at 5.3.7 or 5.3.8 START TO CHANGE SYR (Press START)
- 5.3.4 Install MICRO-VALVE (see section 7.1.1/7.1.2). Touch YES to indicate when valve is installed. INST MICRO VALVE (YES)
- 5.3.5 Install 2.5 mL syringe in left hand drive (see 7.2.1-7.2.3). Touch YES to indicate when syringe is installed. INST A SYR AD 2.5 mL (YES)
- 5.3.6 This syringe is not used so any B-syringe may be in the right hand drive for this application. Touch YES. INST ANY B-SYRINGE (YES)
- 5.3.7 Microprocessor will find the zero position of the syringes when start button on hand unit is pressed. **START** → zero set (Press START)
- 5.3.8 The whole system must be primed with liquid and all air removed, to prevent errors in volume due to compression and expansion of entrapped gas. Put aspiration tubing in reagent and delivery tip in empty beaker. Press start button. Syringe will cycle continuously until start button is pressed again. **START** → prime **N** → CONT (Press START)  
endprime → press START
- 5.3.9 MICROLAB 1000 is now ready to run method no. 20. 1 mL reagent has been aspirated. asp 1.0 mL
- 5.3.10 First aliquot will be dispensed when start button is pressed. disp #1 10 ul (Press START)
- 5.3.11 Second to tenth aliquots will be dispensed automatically and syringe refilled automatically.
- 5.3.12 Cycle will repeat as many times as required. Touch NO to enter stop routine.
- 5.3.13 Programme will terminate when YES is touched with possibility to continue when NO is used. **STOP** OR CONTINUE (YES)
- 5.3.14 MICROLAB 1000 ready to function in any mode. **RUN** OR PROG OR AID
- 5.4 Deleting a stored method (in conjunction with 5.2 and 5.3)** DISPLAY
- 5.4.1 Select AID channel. **RUN** OR **PROG** OR **AID** (YES)

5.4.2 Select DELETE.

**REVIEW** OR **DELETE**  
(NO) (YES)

5.4.3 Either all methods or one method can be deleted.  
To delete this exercise enter method # 20.

DEL ALL METHODS? Y/N  
(NO)  
(2,0,YES)

5.4.4 Method 20 now deleted as can be seen by looking at  
the content of location 20. Select REVIEW.

**RUN** OR **PROG** OR **AID**  
(NO) (NO) (YES)

**REVIEW** OR DELETE  
(YES)

5.4.5 Each time YES is touched, location will advance,  
so to location 20. Touch YES again.

20 NO METHOD  
(YES)

5.4.6 Touch NO twice to exit this channel (see flow chart)

**REVIEW** OR **DELETE**  
(NO) (NO)

5.4.7 MICROLAB 1000 ready to function in any mode.

**RUN** OR PROG OR AID

## SECTION 6 - OPERATING PROCEDURE

### 6.1 Symbols used

| Description                | Symbol             | Mode of operation                            | Programmable range  |
|----------------------------|--------------------|----------------------------------------------|---------------------|
| Method number              | n                  |                                              | 1 - 50              |
| Dilution ratio             | x                  | Diluter                                      | 1 : 1 - 1 : 5000    |
| Number of aliquots         | $x_1 - x_{10}$     | Serial diluter                               | 1 : 1 - 1 : 1000    |
| Diluent volume             | y                  | Repeating Dispenser                          | 1 - 100             |
| Reagent volume             | $v_D$              | Diluter/Double diluter                       | 2.5 $\mu$ L - 25 mL |
|                            | $v_R$              | Dispenser / Repeating Dispenser              | 2.5 $\mu$ L - 25 mL |
|                            | $V_{R1} - V_{R10}$ | Double Diluter (Reagent + Sample + Air gaps) | 0 $\mu$ L - 5 mL    |
| Sample volume              | $v_S$              | Serial Dispenser                             | 2.5 $\mu$ L - 5 mL  |
|                            | $V_{S1} - V_{S10}$ | Diluter                                      | 1 $\mu$ L - 5 mL    |
| Total volume               | $v_T$              | Multiple aspiration                          | 1 $\mu$ L - 5 mL    |
| Wash volume                | $v_W$              | Serial Diluter                               | 25 $\mu$ L - 5 mL   |
| Total reagent volume       | $v_{RT}$           | Serum Distributer/Multiple aspiration        | 0 $\mu$ L - 25 mL   |
| A-Syringe volume           | $v_1$              | Serial Dispenser                             | 2.5 $\mu$ L - 25 mL |
| B-Syringe volume           | $v_2$              |                                              | 250 $\mu$ L - 25 mL |
| Aspiration speed B-syringe | $s_1$              |                                              | 50 $\mu$ L - 5 mL   |
| Dispense speed B-syringe   | $s_2$              |                                              |                     |
| Aspiration speed A-syringe | $s_3$              |                                              |                     |
| Dispense speed A-syringe   | $s_4$              |                                              | 1 - 20 sec.         |

## 6.2 Program as diluter

| Display                                                       | Action                           | Result                                       | Alternative Action | Alternative Result Comment                                      |
|---------------------------------------------------------------|----------------------------------|----------------------------------------------|--------------------|-----------------------------------------------------------------|
| SELF TEST COMPLETED                                           | Switch on power                  | -                                            | Comment            | BRIEFLY DISPLAYED                                               |
| 6.2.1 <b>[RUN]</b> OR PROG OR AID                             | NO                               | Reject RUN channel                           | YES                | RUN channel selected                                            |
| 6.2.2 RUN OR <b>[PROG]</b> OR AID                             | YES                              | Select programme channel                     | NO                 | PROG channel rejected                                           |
| 6.2.3 METHOD #                                                | Enter unoccupied location number | Number displayed                             |                    |                                                                 |
| 6.2.4 METHOD (n)                                              | YES                              | Location referenced for method being entered | NO                 | No effect                                                       |
| 6.2.5 <b>[DIL]</b> DDL DISP DIST                              | YES                              | Select DIL function                          | NO                 | Reject DIL function                                             |
| 6.2.6 <b>[DILUTER]</b> / SERIAL DIL                           | YES                              | Select DIL function                          | NO                 | Select SERIAL DIL function.                                     |
| 6.2.7 AIR GAP 10 $\mu\text{L}$                                | YES                              | Volume entered in memory                     | NO                 | Enter alternative air gap volume ( $v_A$ ) 0 - 30 $\mu\text{L}$ |
| 6.2.8 autom pick up? Y/N                                      | YES                              | Automatic pick up of diluent                 | NO                 | manual pick up of diluent by activating start button.           |
| 6.2.9 SAMPLE $\mu\text{L}$                                    | Enter volume (x)                 |                                              | NO                 | No effect                                                       |
| 6.2.10 SAMPLE ( $v_S$ ) $\mu\text{L}$<br>0-5000 $\mu\text{L}$ | YES                              | Volume entered in memory                     |                    |                                                                 |

| Display                                                                                                                | Action      | Result                                                                   | Alternative Action | Alternative Comment                                                         |
|------------------------------------------------------------------------------------------------------------------------|-------------|--------------------------------------------------------------------------|--------------------|-----------------------------------------------------------------------------|
| 6.2.11 DILUENT <span style="border: 1px solid black;">µL</span> 1+1÷                                                   | NO          | Reject diluent as a vol.                                                 | YES                | Select diluent as volume                                                    |
| 6.2.12 DILUENT <span style="border: 1px solid black;">µL</span> <span style="border: 1px solid black;">1+</span> 1÷    | NO          | Reject diluent as a ratio<br>1 part sample + (x) parts<br>diluent.       | YES                | Select ratio 1 part sample<br>+ (x) parts diluent.                          |
| 6.2.13 DILUENT <span style="border: 1px solid black;">µL</span> 1 + <span style="border: 1px solid black;">1÷</span>   | YES         | Select diluent as ratio<br>of 1 part sample in (x)<br>parts total volume | NO                 | Reject diluent as a ratio<br>of 1 part sample in (x) parts<br>total volume. |
| 6.2.14 DILUTION 1÷                                                                                                     | Enter ratio |                                                                          |                    |                                                                             |
| 6.2.15 DILUTION 1÷ (x)                                                                                                 | YES         | Ratio of 1 part sample in<br>(x) parts total vol. stored                 | NO                 | No effect                                                                   |
| 6.2.16 VERIFY ?                                                                                                        | YES         | Accepts verification                                                     | NO                 | RUN without verification or<br>modification.                                |
| 6.2.17 <span style="border: 1px solid black;">FUNCT</span> VOL SYR SPEED                                               | YES         | idem                                                                     | NO                 | Selects VOL                                                                 |
| 6.2.18 DILUTER                                                                                                         | YES         | idem                                                                     | NO                 | Selects DIL DDL DISP DIST<br>to restart programming.                        |
| 6.2.19 <span style="border: 1px solid black;">FUNCT</span> <span style="border: 1px solid black;">VOL</span> SYR SPEED | YES         | Verify air gap volume                                                    | NO                 | Change air gap volume                                                       |
| 6.2.20 autom. pick up YES                                                                                              | YES         | verifies automatic pick up                                               | NO                 | Selection of autom. or manual pick<br>up.                                   |
| 6.2.21 SAMPLE (v <sub>S</sub> ) <span style="border: 1px solid black;">µL</span>                                       | YES         | Volume entered in memory                                                 | NO                 | Change sample volume                                                        |
| 6.2.22 DILUTION 1-(x)                                                                                                  | YES         | Ratio entered in memory                                                  | NO                 | Change dilution ratio                                                       |
| 6.2.23 <span style="border: 1px solid black;">FUNCT</span> VOL <span style="border: 1px solid black;">SYR</span> SPEED | YES         | Shows computed A-syringe size                                            | NO                 | Selects SPEED channel                                                       |



| Display                                                         | Action             | Result                                    | Alternative Action | Alternative Comment                                                     |
|-----------------------------------------------------------------|--------------------|-------------------------------------------|--------------------|-------------------------------------------------------------------------|
| 6.2.24 A-SYRINGE AD ( $v_1$ ) $\mu$ L, mL                       | YES                | A-Syringe AD value accepted               | NO                 | Shows next suitable syringe size.                                       |
| 6.2.25 B-SYRINGE DAD ( $v_2$ ) $\mu$ L, mL                      | YES                | B-Syringe DAD entered in memory           | NO                 | idem                                                                    |
| 6.2.26 MICRO(MACRO)-VALVE                                       | YES                | Micro valve memorized                     | NO                 | Shows alternative valve                                                 |
| 6.2.27 FUNCT VOL SYR <b>SPEED</b>                               | YES                | Shows microprocessor selected speeds      | NO                 | VERIFY ?                                                                |
| 6.2.28 asp SAMPLE ( $s_1$ ) sec.<br>( $s_1$ = computed time)    | YES                | Accepts computed time to aspirate sample  | NO                 | Change computed time                                                    |
| 6.2.29 disp SAMPLE ( $s_1$ ) sec.<br>( $s_1$ = computed time)   | YES                | Accepts computed time to dispense sample  | NO                 | idem                                                                    |
| 6.2.30 asp DILUENT ( $s_2$ ) sec.<br>( $s_2$ = computed time)   | YES                | Accepts computed time to aspirate diluent | NO                 | idem                                                                    |
| 6.2.31 disp. DILUENT ( $s_3$ ) sec.<br>( $s_3$ = computed time) | YES                | Accepts computed time to dispense diluent | NO                 | idem.                                                                   |
| 6.2.32 VERIFY ?                                                 | NO                 | Selects RUN channel                       | YES                | FUNCT VOL SYR SPEED<br>see 6.2.17                                       |
| 6.2.33 <b>RUN</b> OR PROG OR AID                                | YES                | Selects RUN channel                       | NO                 | Rejects RUN channel                                                     |
| 6.2.34 METHOD # (n)                                             | YES                | Accepts method (n) just programmed        | NO                 | Rejects to run method # (n) and enables entering of another method #.   |
| 6.2.35 (n) DIL ( $v_5$ ) $\mu$ L 1 ÷ (x)                        | YES                | Displays method with volumes              | NO                 | RUN OR PROG OR AID                                                      |
| 6.2.36 START TO CHANGE SYR                                      | Press start button | Lower drives to allow syringe fitting     | Comment            | If syringes already in place are ok, prog continues at 6.2.40 or 6.2.41 |



| Display                                                             | Action                                           | Result                                                                                                | Alternative Action | Alternative Result Comment                                               |
|---------------------------------------------------------------------|--------------------------------------------------|-------------------------------------------------------------------------------------------------------|--------------------|--------------------------------------------------------------------------|
| 6.2.37 INST MICRO(MACRO)-VALVE                                      | Install valve and indicate by YES                | Microprocessor memorizes which valve fitted                                                           | Comment/NO         | Only ever change valve when indicated to do so/NO - continue at 6.2.26   |
| 6.2.38 INST ASYR AD ( $v_1$ ) $\mu\text{L}$                         | Install A-Syringe indicate by YES                | Microprocessor memorizes which A-syringe fitted                                                       | Comment/NO         | Only ever change syringe when indicated to do so/NO - continue at 6.2.24 |
| 6.2.39 INST BSYR DAD ( $v_2$ ) $\mu\text{L}$                        | Install B-Syringe indicate by YES                | Microprocessor memorizes which B-syringe fitted                                                       | Comment/NO         | Only ever change syringe when indicated to do so/NO - continue at 6.2.24 |
| 6.2.40 START $\rightarrow$ zero set                                 | Press start button                               | MICROLAB 1000 finds zero position of syringes                                                         |                    |                                                                          |
| 6.2.41 <b>START</b> $\rightarrow$ prime <b>N</b> $\rightarrow$ cont | Press start button                               | A-syringe cycles continuously until start button is pressed                                           | NO                 | No prime cycle required                                                  |
| 6.2.42 endprime PRESS START                                         | Press start button                               | Prime cycle stops                                                                                     |                    |                                                                          |
| 6.2.43 asp AIR GAP ( $v_A$ ) $\mu\text{L}$                          | automatically                                    | B-Syringe aspirates ( $v_a$ ) $\mu\text{L}$ of air                                                    | YES/NO Comment     | Enters stop routine If no air gap, continue at 6.2.44                    |
| 6.2.44 asp ( $v_D$ ) + ( $v_S$ ) $\mu\text{L}$                      | Put probe tip in sample and press start button   | B-Syringe aspirates sample and A-syringe aspirates ( $v_D$ ) $\mu\text{L}$ diluent from stock diluent | YES/NO             | Enters stop routine                                                      |
| 6.2.45 disp ( $v_D$ ) + ( $v_S$ ) $\mu\text{L}$                     | Put probe tip in receptacle + press start button | Diluted sample is dispensed                                                                           | YES/NO             | Enters stop routine                                                      |
| 6.2.46 As 6.2.43 above                                              | continuous operation                             |                                                                                                       | Comment            | Operation repeats until YES/NO is touched                                |

| Display                          | Action | Result                                       | Alternative Action | Alternative Result Comment |
|----------------------------------|--------|----------------------------------------------|--------------------|----------------------------|
| 6.2.47 As 6.2.44 above           | NO/YES | Stops routine                                | NO                 | Stops routine              |
| 6.2.48 <b>STOP</b> OR CONTINUE   | YES    | DILUTER terminated                           | NO                 | Programme can continue     |
| 6.2.49 <b>RUN</b> OR PROG OR AID | -      | MICROLAB 1000 ready to function in any mode. | -                  | -                          |

### 6.3 Program as double diluter

|        | Display                            | Action                           | Result                                  | Alternative Action | Alternative Result Comment                                       |
|--------|------------------------------------|----------------------------------|-----------------------------------------|--------------------|------------------------------------------------------------------|
|        | SELF TEST COMPLETED                | Switch on power                  | -                                       | Comment            | BRIEFLY DISPLAYED                                                |
| 6.3.1  | <b>RUN</b> OR PROG OR AID          | NO                               | Reject RUN channel                      | YES                | RUN channel selected                                             |
| 6.3.2  | RUN OR <b>PROG</b> OR AID          | YES                              | Select PROG channel                     | NO                 | Reject PROG channel                                              |
| 6.3.3  | Method #                           | enter unoccupied location number | Number displayed                        |                    |                                                                  |
| 6.3.4  | METHOD (n)                         | YES                              | Location referenced for method entering | NO                 | No effect                                                        |
| 6.3.5  | <b>DIL</b> DDL DISP DIST           | NO                               | Reject DIL function                     | YES                | Select DIL function                                              |
| 6.3.6  | DIL <b>DDL</b> DISP DIST           | YES                              | Select DDL function                     | NO                 | Reject DDL function                                              |
| 6.3.7  | AIR GAP 10 $\mu\text{L}$           | YES                              | 1st air gap volume ( $v_A$ ) memorized  | NO                 | Select alternative air gap volume ( $v_A = 0 - 30 \mu\text{L}$ ) |
| 6.3.8  | autom. pick up? Y/N                | YES                              | Automatic pick up of diluent            | NO                 | manual pick up of diluent by pushing the handactuator button     |
| 6.3.9  | REAGENT $\mu\text{L}$              | enter volume                     |                                         |                    |                                                                  |
| 6.3.10 | REAGENT ( $v_R$ ) $\mu\text{L}$    | YES                              | Reagent volume memorized                |                    |                                                                  |
| 6.3.11 | AIR GAP ( $v_{A2}$ ) $\mu\text{L}$ | enter volume 0-30                |                                         | NO                 | No effect                                                        |

| Display                                   | Action       | Result                                          | Alternative Action | Alternative Result Comment                                |
|-------------------------------------------|--------------|-------------------------------------------------|--------------------|-----------------------------------------------------------|
| 6.3.12 AIR GAP ( $v_{A2}$ ) $\mu\text{L}$ | YES          | 2nd air gap vol. ( $v_{A2}$ ) memorized         | NO                 | No effect                                                 |
| 6.3.13 SAMPLE $\mu\text{L}$               | enter volume |                                                 |                    |                                                           |
| 6.3.14 SAMPLE ( $v_S$ ) $\mu\text{L}$     | YES          | Sample vol. memorized                           | NO                 | No effect                                                 |
| 6.3.15 DILUENT $\mu\text{L}$              | enter volume |                                                 |                    |                                                           |
| 6.3.16 DILUENT ( $v_D$ ) $\mu\text{L}$    | YES          | Diluent vol. memorized                          | NO                 | No effect                                                 |
| 6.3.17 VERIFY?                            | YES          | Opportunity to verify or amend parameters       | NO                 | Accept all parameters and select RUN channel              |
| 6.3.18 <b>[FUNCT]</b> VOL SYR SPEED       | YES          | Show function selected                          | NO                 | Rejects verification of function                          |
| 6.3.19 DOUBLE DILUTER                     | YES          | Accepts function                                | NO                 | Rejects function as incorrect. Program continues at 6.3.5 |
| 6.3.20 <b>[VOL]</b> SYR SPEED             | YES          | Displays air gap volume                         | NO                 | Selects SYR                                               |
| 6.3.21 AIR GAP ( $v_{A1}$ ) $\mu\text{L}$ | YES          | Memorizes air gap volume                        | NO                 | Change air gap volume                                     |
| 6.3.22 autom. pick up $\rightarrow$ YES   | YES          | Memorizes autom. pick up of diluent             | NO                 | Option to select manual or automatic pick up.             |
| 6.3.23 Reagent ( $v_R$ ) $\mu\text{L}$    | YES          | Microprocessor memorizes reagent vol. ( $v_R$ ) | NO                 | Change reagent vol. ( $v_R$ )                             |
| 6.3.24 Air gap ( $v_{A2}$ ) $\mu\text{L}$ | YES          | Microprocessor memorizes 2nd air gap vol.       | NO                 | Change 2nd air gap vol. ( $v_{A2}$ )                      |
| 6.3.25 Sample ( $v_S$ ) $\mu\text{L}$     | YES          | Microprocessor memorizes sample vol. ( $v_S$ )  | NO                 | Change sample vol. ( $v_S$ )                              |
| 6.3.26 Diluent ( $v_D$ ) $\mu\text{L}$    | YES          | Microprocessor memorizes diluent vol. ( $v_D$ ) | NO                 | Change diluent vol. ( $v_D$ )                             |

| Display                                                        | Action | Result                                       | Alternative Action | Alternative Result Comment                 |
|----------------------------------------------------------------|--------|----------------------------------------------|--------------------|--------------------------------------------|
| 6.3.27 FUNCT VOL <b>[SYR]</b> SPEED                            | YES    | Shows microprocessor selected syringes.      | NO                 | Rejects SYR channel                        |
| 6.3.28 A-SYRINGE AD ( $v_1$ ) $\mu$ L, mL                      | YES    | Microprocessor memorizes A-syringe           | NO                 | Selects next larger suitable syringe size. |
| 6.3.29 B-SYRINGE DAD ( $v_2$ ) $\mu$ L, mL                     | YES    | Microprocessor memorizes B-syringe.          | NO                 | idem.                                      |
| 6.3.30 MICRO(MACRO)-VALVE                                      | YES    | Display shows suitable valve block           | NO                 | Selects alternative valve                  |
| 6.3.31 FUNCT VOL SYR <b>[SPEED]</b>                            | YES    | Displays computerized speeds                 | NO                 | Selects VERIFY?                            |
| 6.3.32 asp SAMPLE ( $s_1$ ) sec.<br>( $s_1$ = computed speed)  | YES    | Accepts calculated time to aspirate sample   | NO                 | Change calculated time                     |
| 6.3.33 disp SAMPLE ( $s_2$ ) sec<br>( $s_2$ = computed speed)  | YES    | Accepts calculated time to dispense sample   | NO                 | idem.                                      |
| 6.3.34 asp DILUENT ( $s_3$ ) sec<br>( $s_3$ = computed speed)  | YES    | Accepts calculated time to aspirate diluent  | NO                 | idem.                                      |
| 6.3.35 disp DILUENT ( $s_4$ ) sec<br>( $s_4$ = computed speed) | YES    | Accepts calculated time to dispense diluent. | NO                 | idem.                                      |
| 6.3.36 VERIFY ?                                                | NO     | Selects <b>[RUN]</b> channel                 | YES                | Repeats from 6.3.18.                       |
| 6.3.37 RUN OR PROG OR AID                                      | YES    | Selects channel to run program.              | NO                 | Rejects last method.                       |
| 6.3.38 LAST METHOD # (n)                                       | YES    | Displays method just programmed.             | NO                 | Rejects last method.                       |
| 6.3.39 (n) DDL ( $v_R$ )( $v_S$ ) + ( $v_D$ )                  | YES    | Displays method with volumes                 | NO                 | RUN OR PROG OR AID.                        |

| Display                                                             | Action                                        | Result                                                       | Alternative Action | Alternative Result Comment                                                        |
|---------------------------------------------------------------------|-----------------------------------------------|--------------------------------------------------------------|--------------------|-----------------------------------------------------------------------------------|
| 6.3.40 START TO CHANGE SYR                                          | Press start button                            | Lower drives to allow fitting of syringes                    | Comment            | If syringes fitted are ok the program will continue at 6.3.44 or 6.3.45           |
| 6.3.41 INST MICRO(MACRO)-VALVE                                      | Install valve indicate by YES                 | Computer memorizes which valve is fitted                     | Comment/NO         | Only ever change valve when instructed by the computer/NO -continue at 6.3.30.    |
| 6.3.42 INST ASYR AD ( $v_1$ ) $\mu\text{L}$                         | Install A-syringe indicate by YES             | Computer memorizes A-syringe ( $v_1$ ) $\mu\text{L}$ fitted  | Comment/NO         | Only ever change syringes when instructed by the computer/NO -continue at 6.3.28. |
| 6.3.43 INST BSYR DAD ( $v_2$ ) $\mu\text{L}$                        | Install B-syringe indicate by YES             | Computer memorizes B-syringe ( $v_2$ ) $\mu\text{L}$ fitted  | Comment/NO         | idem.                                                                             |
| 6.3.44 <b>START</b> $\rightarrow$ zero set                          | Press start button                            | MICROLAB 1000 finds zero position of syringes                | Comment            | If zero already set, continue 6.3.45.                                             |
| 6.3.45 <b>START</b> $\rightarrow$ prime <b>N</b> $\rightarrow$ cont | Press start button                            | A-Syringe cycles continuously until start button is pressed. | Comment            | Prime liquid system until no air bubbles are visible.                             |
| 6.3.46 endprime - PRESS START                                       | Press start button                            | Prime cycle stops.                                           |                    |                                                                                   |
| 6.3.47 asp AIR GAP ( $v_{A1}$ ) $\mu\text{L}$                       | asp automatically                             | B-Syringe aspirates ( $v_{A1}$ ) $\mu\text{L}$ air gap       |                    |                                                                                   |
| 6.3.48 asp ( $v_D$ ) + ( $v_R$ ) $\mu\text{L}$                      | Press start button place probe tip in reagent | B-Syringe aspirates reagent and A-syringe aspirates diluent. | NO                 | STOP OR CONTINUE                                                                  |

| Display                                               | Action                                                      | Result                                         | Alternative Action | Alternative Result Comment |
|-------------------------------------------------------|-------------------------------------------------------------|------------------------------------------------|--------------------|----------------------------|
| 6.3.49 asp AIR GAP ( $v_{A2}$ ) $\mu\text{L}$         | Press start button<br>remove probe tip<br>from reagent      | B-Syringe aspirates 2nd<br>air gap             | NO                 | idem.                      |
| 6.3.50 asp ( $v_S$ ) $\mu\text{L}$                    | Press start button<br>place probe tip in<br>sample          | B-Syringe aspirates sample                     | NO                 | idem.                      |
| 6.3.51 disp ( $v_D$ ) + ( $v_R + v_S$ ) $\mu\text{L}$ | Put delivery tip<br>into receptacle +<br>press start button | Diluted sample delivered                       | NO                 | idem.                      |
| 6.3.52 As 6.3.49 above                                |                                                             | Cycle repeats until NO/YES<br>is touched.      | NO                 | idem.                      |
| 6.3.53 As 6.3.50 above                                | NO/YES                                                      | Enters stop routine                            |                    |                            |
| 6.3.54 <b>STOP</b> OR CONTINUE                        | YES                                                         | Terminates diluter                             | NO                 | Program continues          |
| 6.3.55 <b>RUN</b> OR PROG OR AID                      | -                                                           | MICROLAB 1000 ready to<br>function in any mode | -                  | -                          |



#### 6.4 Program as single or multiple dispenser

| Display                           | Action                           | Result                                       | Alternative Action | Alternative Result Comment                                                          |
|-----------------------------------|----------------------------------|----------------------------------------------|--------------------|-------------------------------------------------------------------------------------|
| SELF TEST COMPLETED               | Switch on power                  | -                                            | -                  | BRIEFLY DISPLAYED                                                                   |
| 6.4.1 <b>RUN</b> OR PROG OR AID   | NO                               | Reject RUN channel                           | YES                | Select RUN channel                                                                  |
| 6.4.2 RUN OR <b>PROG</b> OR AID   | YES                              | Select PROG channel                          | NO                 | Reject PROG channel                                                                 |
| 6.4.3 METHOD #                    | enter unoccupied location number | Number displayed                             |                    |                                                                                     |
| 6.4.4 METHOD # (n)                | YES                              | Location referenced for new program.         | NO                 | No effect                                                                           |
| 6.4.5 <b>DIL</b> DDL DISP DIST    | NO                               | Reject DIL function                          | YES                | Select DIL function                                                                 |
| 6.4.6 DIL <b>DDL</b> DISP DIST    | NO                               | Reject DDL function                          | YES                | Select DDL function                                                                 |
| 6.4.7 DIL DDL <b>DISP</b> DIST    | YES                              | Select DISP function                         | NO                 | Reject DISP function                                                                |
| 6.4.8 <b>REP.DISP</b> SERIAL DISP | YES                              | Select REP DISP function                     | NO                 | Reject REP DISP function                                                            |
| 6.4.9 asp.autom.Σ vol? Y/N        | YES                              | Select automatic aspiration                  | NO                 | Reject automatic aspiration, but manual operation by pressing start button          |
| 6.4.10 disp. automatic? Y/N       | NO                               | Activate dispensing by pressing start button | YES                | Automatic dispensing of each aliquot, with selectable time delay of 0.1 to 9.9 sec. |

| Display                                          | Action                | Result                                                 | Alternative Action | Alternative Result Comment          |
|--------------------------------------------------|-----------------------|--------------------------------------------------------|--------------------|-------------------------------------|
| 6.4.11 REAGENT $\mu\text{L}$                     | enter volume          | Volume is displayed in $\mu\text{L}$                   |                    |                                     |
| 6.4.12 REAGENT ( $v_R$ ) $\mu\text{L}$           | YES                   | Volume stored in memory                                | NO                 | No effect                           |
| 6.4.13 # OF ALIQUOTS                             | enter no. of aliquots | Shows number of aliquots                               | NO                 | No effect                           |
| 6.4.14 # OF ALIQUOTS (y)                         | YES                   | Computer calculates volume of total number of aliquots | NO                 | No effect                           |
| 6.4.15 VERIFY?                                   | YES                   | Opportunity to verify or amend program                 | NO                 | Accept all parameters               |
| 6.4.16 <b>FUNCT</b> VOL SYR SPEED                | YES                   | Displays selected method of operation                  | NO                 | Selects VOL                         |
| 6.4.17 REPEATING DISPENSER                       | YES                   | Memorizes selected mode of operation                   | NO                 | Selects SERIAL DISPENSER            |
| 6.4.18 asp. automatic $\Sigma$ VOL               | YES                   | Memorizes automatic aspiration of total reagent volume | NO                 | Manual operation                    |
| 6.4.18.1 manual START dispense                   | YES                   |                                                        | NO                 | automatic dispense                  |
| 6.4.18.2 <b>FUNCT</b> VOL SYR SPEED              | YES                   | Displays volumes                                       | NO                 | Reject VOL channel                  |
| 6.4.19 REAGENT ( $v_R$ ) $\mu\text{L}$           | YES                   | Displays reagent volume                                | NO                 | Change reagent volume               |
| 6.4.20 # OF ALIQUOTS (y)                         | YES                   | Memorizes (y) # of aliquots                            | NO                 | Change # of aliquots                |
| 6.4.21 <b>FUNCT</b> VOL <b>SYR</b> SPEED         | YES                   | Displays computed syringe sizes                        | NO                 | Rejects SYR channel                 |
| 6.4.22 A-SYRINGE AD ( $v_1$ ) $\mu\text{L}$ , mL | YES                   | Memorizes calculated syringe size                      | NO                 | Displays next suitable syringe size |

| Display                                              | Action                            | Result                                                                | Alternative Action | Alternative Result Comment                                                            |
|------------------------------------------------------|-----------------------------------|-----------------------------------------------------------------------|--------------------|---------------------------------------------------------------------------------------|
| 6.4.23 MICRO(MACRO)-VALVE                            | YES                               | Micro-valve accepted                                                  | NO                 | Shows alternative valve                                                               |
| 6.4.24 FUNCT VOL SYR <b>SPEED</b>                    | YES                               | Shows microprocessor selected speeds                                  | NO                 | VERIFY?                                                                               |
| 6.4.25 asp REAGENT (s3) sec<br>(s3 = computed time)  | YES                               | Accept calculated time to aspirate reagent                            | NO                 | Change calculated time                                                                |
| 6.4.26 disp REAGENT (s4) sec<br>(s4 = computed time) | YES                               | Accept calculated time to dispense reagent                            | NO                 | idem                                                                                  |
| 6.4.27 VERIFY?                                       | NO                                | Selects RUN channel                                                   | YES                | Verification required, repetition from step 6.4.16                                    |
| 6.4.28 <b>RUN</b> OR PROG OR AID                     | YES                               | RUN channel selected                                                  | NO                 | Reject RUN channel                                                                    |
| 6.4.29 METHOD # (n)                                  | YES                               | Displays methods (n)                                                  |                    |                                                                                       |
| 6.4.30 (n) RDISP man (y) x (v <sub>R</sub> ) µL      | YES                               | Accepts method (n) just programmed and displays description of method | NO                 | Select other method number                                                            |
| 6.4.31 START TO CHANGE SYR                           | Press start button                |                                                                       | Comment            | If correct syringes are already installed, program continues at 6.4.35 START zero set |
| 6.4.32 INST MICRO(MACRO)-VALVE                       | Install valve and indicate by YES | Microprocessor memorizes which valve fitted                           | Comment NO         | Only ever change valve when indicated to do so/NO - continue at 6.4.23                |
| 6.4.33 INST ASYR AD                                  | Install A-syringe indicate by YES | Microprocessor memorizes which A-syringe fitted                       | Comment NO         | Only ever change syringe when indicated to do so/NO - continue at 6.4.22              |

| Display                                     | Action                       | Result                                                            | Alternative Action | Alternative Result Comment                      |
|---------------------------------------------|------------------------------|-------------------------------------------------------------------|--------------------|-------------------------------------------------|
| 6.4.34 INST ANY B-SYRINGE                   | YES                          | Size of syringe has no impact on results                          | NO                 | Continue to START zero set                      |
| 6.4.35 <b>START</b> → zero set              | Press start button           | MICROLAB 1000 finds zero position of syringes                     | NO                 | Press start button                              |
| 6.4.36 <b>START</b> → prime <b>N</b> → cont | Press start button           | A-Syringe cycles continuously until start button is pressed again | NO                 | Program continues without priming               |
| 6.4.37 endprime → PRESS START               | Press start button           | Prime cycle stops.                                                |                    |                                                 |
| 6.4.38 asp ( $v_R$ ) $\mu$ L                | Press start button           | Aspiration of total volume from stock solution                    |                    |                                                 |
| 6.4.39 disp #1 ( $v_R$ ) $\mu$ L            | Press start button           | First aliquot ( $v$ ) is dispensed                                | NO                 | STOP OR CONTINUE                                |
| 6.4.40 disp # 2 ... (y) x ( $v_R$ ) $\mu$ L | Press start button           | 2nd and subsequent aliquots are dispensed                         | NO                 | idem                                            |
| 6.4.41                                      | Press start button<br>NO/YES | Program is interrupted                                            |                    |                                                 |
| 6.4.42 <b>STOP</b> OR CONTINUE              | YES                          | Program terminated                                                | NO                 | Program continues at position where interrupted |
| 6.4.43 <b>RUN</b> OR PROG OR AID            | -                            | MICROLAB 1000 ready to function in any mode                       |                    |                                                 |

## 6.5 Program as a serum distributor

| Display                           | Action                           | Result                              | Alternative Action | Alternative Result Comment                             |
|-----------------------------------|----------------------------------|-------------------------------------|--------------------|--------------------------------------------------------|
| SELF TEST COMPLETED               | Switch on power                  | -                                   | Comment            | BRIEFLY DISPLAYED                                      |
| 6.5.1 <b>RUN</b> OR PROG OR AID   | NO                               | Reject RUN channel                  | YES                | Select RUN channel                                     |
| 6.5.2 RUN OR <b>PROG</b> OR AID   | YES                              | Program channel selected            | NO                 | Reject PROG channel                                    |
| 6.5.3 Method #                    | Enter unoccupied location number | Number displayed                    |                    |                                                        |
| 6.5.4 METHOD # (n)                | YES                              | Location referenced for new program | NO                 | No effect                                              |
| 6.5.5 <b>DIL</b> DDL DISP DISTR   | NO                               | Reject DIL function                 | YES                | Select DIL function                                    |
| 6.5.6 DIL <b>DDL</b> DISP DIST    | NO                               | Reject DDL function                 | YES                | Select DDL function                                    |
| 6.5.7 DIL DDL <b>DISP</b> DIST    | NO                               | Reject DISP function                | YES                | Select DISP function                                   |
| 6.5.8 DIL DDL DISP <b>DIST</b>    | YES                              | Select DIST function                | NO                 | Reject DIST function                                   |
| 6.5.9 <b>SERUM DIST</b> MULTI ASP | YES                              | Select SERUM DIST function          | NO                 | Selects MULTI ASP                                      |
| 6.5.10 AIR GAP 30 µL              | YES                              | Air gap volume confirmed            | NO                 | Automatic pick up? Y/N                                 |
| 6.5.11 automatic pick up? Y/N     | YES                              | Automatic pick up of wash volume    | NO                 | Manual pick up of wash volume by pressing start button |

| Display                                  | Action                                                     | Result                                        | Alternative Action | Alternative Result Comment                             |
|------------------------------------------|------------------------------------------------------------|-----------------------------------------------|--------------------|--------------------------------------------------------|
| 6.5.12 SAMPLE # (v <sub>S1</sub> ) µL    | enter volume                                               | Volume displayed                              | NO                 | Terminates sample volume entry                         |
| 6.5.13 SAMPLE # (v <sub>S1</sub> ) µL    | YES                                                        | Volume memorized                              | NO                 | No effect                                              |
| 6.5.14 SAMPLE # (v <sub>S2</sub> ) µL    | enter volume                                               | Volume displayed                              | NO                 | Terminates sample volume entry.<br>Selects wash volume |
| .                                        |                                                            |                                               |                    |                                                        |
| .                                        |                                                            |                                               |                    |                                                        |
| .                                        |                                                            |                                               |                    |                                                        |
| 6.5.15 SAMPLE # (v <sub>S10</sub> ) µL   | enter volume                                               | Volume displayed                              | NO                 | Terminates sample volume entry.<br>Selects wash volume |
| 6.5.16 Sample # (v <sub>S10</sub> ) µL   | YES                                                        | Volume memorized                              | NO                 | idem                                                   |
| 6.5.17 WASH VOLUME µL                    | enter volume<br>(recommend volume:<br>10xtotal serum vol.) | Volume displayed                              | NO                 | No effect                                              |
| 6.5.18 VERIFY?                           | YES                                                        | Opportunity to verify or<br>amend program     | NO                 | Selects RUN channel                                    |
| 6.5.19 <b>FUNCT</b> VOL SYR SPEED        | YES                                                        | Shows mode of operation                       | NO                 | Selects VOL channel                                    |
| 6.5.20 SERUM DISTRIBUTOR                 | YES                                                        | Memorizes function programmed                 | NO                 | Selects DIL DDL DISP DIST                              |
| 6.5.21 <b>FUNCT</b> <b>VOL</b> SYR SPEED | YES                                                        | Displays programmed volumes                   | NO                 | Selects SYR channel                                    |
| 6.5.22 AIR GAP 30 µL                     | YES                                                        | Displays mode of pick up of<br>wash solution  |                    |                                                        |
| 6.5.23 autom.pick up → YES               | YES                                                        | Memorizes mode of pick up<br>of wash solution | NO                 | Automatic pick ? Y/N                                   |

| Display                                                      | Action | Result                                      | Alternative Action | Alternative Result Comment       |
|--------------------------------------------------------------|--------|---------------------------------------------|--------------------|----------------------------------|
| 6.5.24 SAMPLE # 1 ( $v_{S1}$ ) $\mu$ L                       | YES    | Memorizes 1st sample volume                 | NO                 | Change 1st sample volume         |
| 6.5.25 SAMPLE # 2 ( $v_{S2}$ ) $\mu$ L                       | YES    | Memorizes 2nd sample volume                 | NO                 | Change 2nd sample volume         |
| 6.5.26 SAMPLE # 3 ( $v_{S3}$ ) $\mu$ L                       | YES    | Memorizes 3rd sample volume                 | NO                 | Change 3rd sample volume         |
| 6.5.27 SAMPLE # 10 ( $v_{S10}$ ) $\mu$ L                     | YES    | Memorizes 10th sample volume                | NO                 | Change 10th sample volume        |
| 6.5.28 WASH VOLUME ( $v_w$ ) $\mu$ L                         | YES    | Memorizes wash volume                       | NO                 | Change wash volume               |
| 6.5.29 FUNCT VOL <b>SYR</b> SPEED                            | YES    | Shows microprocessor selected syringe sizes | NO                 | Rejects syringe verification     |
| 6.5.30 A-SYRINGE AD ( $v_1$ ) mL                             | YES    | A-Syringe value accepted                    | NO                 | Shows next suitable syringe size |
| 6.5.31 B-SYRINGE DAD ( $v_2$ ) $\mu$ L, mL                   | YES    | B-Syringe value accepted                    | NO                 | Shows next suitable syringe size |
| 6.5.32 MICRO(MACRO)-VALVE                                    | YES    | Valve type accepted                         | NO                 | Shows alternative valve type.    |
| 6.5.33 FUNCT VOL SYR <b>SPEED</b>                            | YES    | Shows computed speeds                       | NO                 | VERIFY ?                         |
| 6.5.34 asp SAMPLE ( $s_3$ ) sec<br>( $s_3$ = computed time)  | YES    | Accept computed time to aspirate sample     | NO                 | Change computed time.            |
| 6.5.35 disp SAMPLE ( $s_4$ ) sec<br>( $s_4$ = computed time) | YES    | Accept computed time dispense sample        | NO                 | idem                             |



| Display                                                             | Action                            | Result                                          | Alternative Action | Alternative Result Comment                                                       |
|---------------------------------------------------------------------|-----------------------------------|-------------------------------------------------|--------------------|----------------------------------------------------------------------------------|
| 6.5.36 asp WASH ( $s_1$ ) sec ( $s_1$ = computed time)              | YES                               | Accept computed time to aspirate wash solution  | NO                 | idem                                                                             |
| 6.5.37 disp WASH ( $s_2$ ) sec ( $s_2$ = computed time)             | YES                               | Accept computed time to dispense wash solution  | NO                 | idem                                                                             |
| 6.5.38 VERIFY?                                                      | NO                                | All parameters accepted                         | NO                 | Verification required                                                            |
| 6.5.39 <b>RUN</b> OR PROG OR AID                                    | YES                               | Select RUN channel to execute method programmed | NO                 | Selects PROG channel                                                             |
| 6.5.40 METHOD # (n)                                                 | YES                               | Accepts method (n) just programmed              | NO                 | Select another method                                                            |
| 6.5.41 (n) DIST $\Sigma$ (v waste + $v_{S1}$ to $v_{S10}$ ) + $v_W$ |                                   | Displays method and volumes                     | NO                 | RUN OR PROG OR AID                                                               |
| 6.5.42 START TO CHANGE SYR                                          | Press start button                |                                                 | Comment            | If correct syringes are already installed, program continues at 6.5.45 or 6.5.46 |
| 6.5.43 INST MICRO(MACRO)-VALVE                                      | Install valve indicate by YES     | Computer memorizes which valve is fitted        | Comment/NO         | Only ever change valve when indicated to do so/NO - continue at 6.5.31           |
| 6.5.44 INST ASYR AD                                                 | Install A-syringe indicate by YES | Computer memorizes which A-syringe fitted       | Comment/NO         | Only ever change syringes when indicated to do so/NO - continue at 6.5.29        |
| 6.5.45 INST BSYR DAD                                                | Install B-syringe indicate by YES | Computer memorizes which B-syringe fitted       | Comment/NO         | Only ever change syringes when indicated to do so/NO - continue at 6.5.30        |
| 6.5.46 <b>START</b> → zero set                                      | Press start button                | MICROLAB 1000 finds zero position of syringes   | NO                 | Press start button                                                               |

| Display                                                         | Action                                            | Result                                                | Alternative Action | Alternative Result Comment |
|-----------------------------------------------------------------|---------------------------------------------------|-------------------------------------------------------|--------------------|----------------------------|
| 6.5.47 <b>START</b> →prime <b>N</b> →cont                       | NO                                                | Priming of system not required                        | YES                | priming required           |
| 6.5.48 asp AIR GAP 30 $\mu$ L                                   | Press start button                                | B-syringe aspirates air gap                           | NO                 | STOP OR CONTINUE           |
| 6.5.49 asp WASH ( $v_w$ ) + SAMPLE<br>( $v_{S1}$ to $v_{S10}$ ) | Press start button<br>Place tip in sample         | Aspiration of total wash and sample volume            | NO                 | idem                       |
| 6.5.50 disp # 1 ( $v_{S1}$ ) $\mu$ L                            | Press start button<br>Place tip in 1st receptacle | 1st sample delivered                                  | YES/NO             | idem                       |
| 6.5.51 disp # 2 ( $v_{S2}$ ) $\mu$ L                            | idem<br>Place tip in 2nd receptacle               | 2nd sample delivered                                  | YES/NO             | idem                       |
| 6.5.52 disp # 3 ( $v_{S3}$ ) $\mu$ L                            | idem<br>Place tip in 3rd receptacle               | 3rd sample delivered                                  | YES/NO             | idem                       |
| 6.5.53 disp # 10 ( $v_{S10}$ ) $\mu$ L                          | idem<br>Place tip in 10th receptacle              | 10th sample delivered                                 | YES/NO             | idem                       |
| 6.5.54 disp WASTE + wash                                        | idem<br>Place tip in receptacle                   | Extra sample and tip washed with ( $v$ ) $\mu$ L sol. | YES/NO             | idem                       |
| 6.5.55 As 6.5.47 above                                          | idem                                              | Cycle repeats                                         | YES/NO             | idem                       |

| Display                          | Action | Result                                      | Alternative Action | Alternative Result Comment          |
|----------------------------------|--------|---------------------------------------------|--------------------|-------------------------------------|
| 6.5.56 As 6.6.48 above           | NO     | Enters stop routine                         | START              | Cycle continues                     |
| 6.5.57 <b>STOP</b> OR CONTINUE   | YES    | Terminates distributor                      | NO                 | Program continues where interrupted |
| 6.5.58 <b>RUN</b> OR PROG OR AID | -      | MICROLAB 1000 ready to function in any mode | -                  | -                                   |

6.5.59 MICROLAB 1000 ready to function in any mode

# 6.6. Program as serial diluter

| Display                                  | Action                           | Result                       | Alternative Action | Alternative Result Comment                           |
|------------------------------------------|----------------------------------|------------------------------|--------------------|------------------------------------------------------|
| SELF TEST COMPLETED                      | Switch power on                  | Microprocessor checks system | Comment            | BRIEFLY DISPLAYED                                    |
| 6.6.1 <span>[RUN]</span> OR PROG OR AID  | NO                               | Rejects RUN channel          | YES                | RUN channel selected                                 |
| 6.6.2 RUN OR <span>[PROG]</span> OR AID  | YES                              | Selects PROG channel         | NO                 | PROG channel rejected                                |
| 6.6.3 METHOD #                           | enter unoccupied location number | Number displayed             | NO                 | No effect                                            |
| 6.6.4 <span>[DIL]</span> DDL DISP DIST   | YES                              | Selects DIL function         | NO                 | Rejects DIL function                                 |
| 6.6.5 <span>[DILUTER]</span> /SERIAL DIL | NO                               | Rejects DIL function         | YES                | Selects DIL function                                 |
| 6.6.6 DILUTER/ <span>[SERIAL DIL]</span> | YES                              | Selects SERIAL DIL function  | NO                 | Rejects SERIAL DIL function                          |
| 6.6.7 AIR GAP 10 $\mu$ L                 | YES                              | Memorizes air gap volume     | NO                 | Change air gap volume                                |
| 6.6.8 autom pick up? Y/N                 | YES                              | Selects automatic pick up    | NO                 | Manual operation, by pressing start button           |
| 6.6.9 TOTAL VOLUME ( $v_T$ ) $\mu$ L     | enter volume of sample + diluent | Volume displayed             | NO                 | No effect                                            |
| 6.6.10 DIL # 1 1: $\div$                 | enter dilution ratio (lowest)    | Lowest dilution ratio        | NO                 | idem<br>VERIFY? : terminates entering dilution ratio |
| 6.6.11 DIL # 1 1: ( $x_1$ )              | YES                              | 1st dilution ratio memorized | NO                 | idem                                                 |

| Display                                | Action               | Result                        | Alternative Action | Alternative Result Comment                                                                                                                                                                    |
|----------------------------------------|----------------------|-------------------------------|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.6.12 DIL # 2 1 ÷                     | enter dilution ratio | 2nd dilution ratio displayed  | NO                 | idem                                                                                                                                                                                          |
| 6.6.13 DIL # 2 1 ÷ (x <sub>2</sub> )   | YES                  | 2nd dilution ratio memorized  | NO                 | idem                                                                                                                                                                                          |
| 6.6.14 DIL # 3 1 ÷                     | enter dilution ratio | 3rd dilution ratio displayed  | NO                 | idem                                                                                                                                                                                          |
| 6.6.15 DIL # 3 1 ÷ (x <sub>3</sub> )   | YES                  | 3rd dilution ratio memorized  | NO                 | idem                                                                                                                                                                                          |
| 6.6.16 DIL # 4 1 ÷                     | enter dilution ratio | 4th dilution ratio displayed  | NO                 | idem                                                                                                                                                                                          |
| 6.6.17 DIL # 4 1 ÷ (x <sub>4</sub> )   | YES                  | 4th dilution ratio memorized  | NO                 | idem                                                                                                                                                                                          |
| 6.6.18 DIL # 10 1 ÷                    | enter dilution ratio | 10th dilution ratio displayed | NO                 | idem                                                                                                                                                                                          |
| 6.6.19 DIL # 10 1 ÷ (x <sub>10</sub> ) | YES                  | 10th dilution ratio memorized | NO                 | idem<br>IMPORTANT<br>Depending on total volume of sample + diluent and dilution ratios selected, certain ratios cannot be executed due to the resolution of the syringe drive. Display shows: |

| Display | Action                             | Result | Alternative Action | Alternative Result Comment                                                       |
|---------|------------------------------------|--------|--------------------|----------------------------------------------------------------------------------|
| 6.6.20  | VERIFY ?                           | YES    | NO                 | DILUTION IMPOSSIBLE<br>NEXT DIL POSSIBLE<br>and closest ratio will be displayed. |
| 6.6.21  | <b>FUNCT</b> VOL SYR SPEED         | YES    | NO                 | Accept with YES and continue, or<br>reject with NO and try different<br>ratio.   |
| 6.6.22  | SERIAL DILUTER                     | YES    | NO                 | RUN OR PROG OR AID                                                               |
| 6.6.23  | <b>FUNCT</b> <b>VOL</b> SYR SPEED  | YES    | NO                 | Selects VOL                                                                      |
| 6.6.24  | AIR GAP ( $v_A$ ) $\mu$ L          | YES    | NO                 | <b>DIL</b> DDL DISP DIST                                                         |
| 6.6.25  | autom. pick up $\rightarrow$ YES   | YES    | NO                 | Selects SYR                                                                      |
| 6.6.26  | TOTAL VOLUME ( $v_T$ ) $\mu$ L, mL | YES    | NO                 | Change air gap volume                                                            |
| 6.6.27  | DIL # 1 1: ( $x_1$ )               | YES    | NO                 | Selects autom pick up? Y/N                                                       |
| 6.6.28  | DIL # 2 1: ( $x_2$ )               | YES    | NO                 | Change volume of diluent<br>and sample                                           |
| 6.6.29  | DIL # 3 1: ( $x_3$ )               | YES    | NO                 | Change dilution ratio                                                            |
|         |                                    |        | NO                 | Change dilution ratio                                                            |
|         |                                    |        | NO                 | Change dilution ratio                                                            |

| Display                                       | Action | Result                                                    | Alternative Action | Alternative Result Comment                                                                       |
|-----------------------------------------------|--------|-----------------------------------------------------------|--------------------|--------------------------------------------------------------------------------------------------|
| 6.6.30 DIL # 10 1: (x <sub>10</sub> )         | YES    | Displays and memorizes 10th and last dilution ratio       | NO                 | Change dilution ratio                                                                            |
| 6.6.31 FUNCT VOL <b>SYR</b> SPEED             | YES    | Displays calculated syringe sizes                         | NO                 | Selects SPEED                                                                                    |
| 6.6.32 A-SYRINGE AD (v <sub>1</sub> ) μL, mL  | YES    | Install requested syringe and confirm by pressing YES     | NO                 | IMPORTANT: no effect, calculated syringe <u>must</u> be installed                                |
| 6.6.33 B-SYRINGE DAD (v <sub>2</sub> ) μL, mL | YES    | Install requested syringe and confirm by pressing YES     | NO                 | idem                                                                                             |
| 6.6.34 MICRO(MACRO)-VALVE                     | YES    | Install requested valve block and confirm by pressing YES | NO                 | Alternative valve block may be installed, as long as no syringes larger than 5 mL are installed. |
| 6.6.35 FUNCT VOL SYR <b>SPEED</b>             | YES    | Computed speeds are displayed                             | NO                 | VERIFY ?                                                                                         |
| 6.6.36 asp SAMPLE (s <sub>3</sub> ) sec       | YES    | Memorizes computed aspiration speed                       | NO                 | Change speed                                                                                     |
| 6.6.37 disp SAMPLE (s <sub>4</sub> ) sec      | YES    | Memorizes computed dispense speed                         | NO                 | Change speed                                                                                     |
| 6.6.38 asp DILUENT (s <sub>1</sub> ) sec      | YES    | Memorizes computed diluent aspiration speed               | NO                 | Change speed                                                                                     |
| 6.6.39 disp DILUENT (s <sub>2</sub> ) sec     | YES    | Memorizes computed diluent dispense speed                 | NO                 | Change speed                                                                                     |



| Display                                                             | Action             | Result                                                                           | Alternative Action | Alternative Result Comment                                                                                         |
|---------------------------------------------------------------------|--------------------|----------------------------------------------------------------------------------|--------------------|--------------------------------------------------------------------------------------------------------------------|
| 6.6.40 VERIFY ?                                                     | NO                 | Selects RUN channel                                                              | YES                | Verifies again from step 6.6.20                                                                                    |
| 6.6.41 <b>RUN</b> OR PROG OR AID                                    | YES                | Selects method # just programmed                                                 | NO                 | RUN OR PROG OR AID                                                                                                 |
| 6.6.42 METHOD # (n)                                                 | YES                | Location referenced and method parameters displayed                              | NO                 | No effect                                                                                                          |
| 6.6.43 (n) SDIL ( $v_T$ ) $\mu$ L 1: (x)                            | YES                | Execution of method confirmed                                                    | NO                 | RUN OR PROG OR AID                                                                                                 |
| 6.6.44 START TO CHANGE SYR                                          | press start button | Lower drives to allow syringe fitting                                            | NO/Comment         | PRESS START<br>If correct syringes are already installed, program continues at <b>START</b> $\rightarrow$ zero set |
| 6.6.45 INST MICRO(MACRO)-VALVE                                      | YES                | Install requested valve block and confirm by YES                                 | NO                 | No effect - requested valve <u>mandatory</u>                                                                       |
| 6.6.46 INST ASYR AD ( $v_1$ ) $\mu$ L, mL                           | YES                | Install requested syringe and confirm by YES                                     | NO                 | No effect - requested syringe must be <u>installed</u>                                                             |
| 6.6.47 INST BSYR DAD ( $v_2$ ) $\mu$ L                              | YES                | Install requested syringe and confirm by YES                                     | NO                 | No effect - requested syringe must be <u>installed</u>                                                             |
| 6.6.48 <b>START</b> zero set                                        | press start button | MICROLAB 1000 finds zero position of syringes                                    | NO                 | PRESS START                                                                                                        |
| 6.6.49 <b>START</b> $\rightarrow$ prime <b>N</b> $\rightarrow$ cont | press start button | A-Syringe cycles continuously until cycle is terminated by pressing start button | NO                 | asp # 1 ( $v_{D1}$ ) + ( $v_{S1}$ )                                                                                |
| 6.6.50 endprime $\rightarrow$ PRESS START                           | press start button | Prime cycle is stopped and volumes of 1st dilution ratio are displayed           | NO                 | <b>STOP</b> OR CONTINUE                                                                                            |

| Display                                     | Action                                                   | Result                                        | Alternative Action | Alternative Result Comment                          |
|---------------------------------------------|----------------------------------------------------------|-----------------------------------------------|--------------------|-----------------------------------------------------|
| 6.6.51 asp # 1 ( $v_{D1}$ ) + ( $v_{S1}$ )  | place probe tip into first sample and press start button | Diluent and sample are aspirated              | NO                 | idem                                                |
| 6.6.52 disp # 1 ( $v_{D1}$ ) + ( $v_{S1}$ ) | press start button to dispense in receptacle             | Diluent and sample are dispensed              | NO                 | idem                                                |
| 6.6.53                                      |                                                          | Dilution series continued until NO is touched | NO                 | Dilution series continue after YES has been touched |
| 6.6.54 <b>STOP</b> OR CONTINUE              | YES                                                      | Enters stop routine                           | NO                 |                                                     |
| 6.6.55 <b>RUN</b> OR PROG OR AID            | -                                                        | MICROLAB 1000 ready to function in any mode   | NO                 |                                                     |

## 6.7 Program as Serial Dispenser

| Display                               | Action                           | Result                                       | Alternative Action | Alternative Result Comment               |
|---------------------------------------|----------------------------------|----------------------------------------------|--------------------|------------------------------------------|
| SELF TEST COMPLETED                   | Switch on power                  | -                                            | Comment            | BRIEFLY DISPLAYED                        |
| 6.7.1 <b>RUN</b> OR PROG OR AID       | NO                               | Reject RUN channel                           | YES                | RUN channel selected                     |
| 6.7.2 RUN OR <b>PROG</b> OR AID       | YES                              | Select PROG channel                          | NO                 | PROG channel rejected                    |
| 6.7.3 Method #                        | enter unoccupied location number | Number displayed                             | -                  | -                                        |
| 6.7.4 Method (n)                      | YES                              | Location referenced for method being entered | NO                 | No effect                                |
| 6.7.5 <b>DIL</b> DDL DISP DIST        | NO                               | Reject DIL function                          | YES                | Select DIL function                      |
| 6.7.6 DIL <b>DDL</b> DISP DIST        | NO                               | Reject DDL function                          | YES                | Select DDL function                      |
| 6.7.7 DIL DDL <b>DISP</b> DIST        | YES                              | Select DISP function                         | YES                | Reject DISP function                     |
| 6.7.8 <b>REP.DISP</b> SERIAL DISP     | NO                               | Reject REP.DISP function                     | YES                | Select REP.DISP function                 |
| 6.7.9 REP.DISP <b>SERIAL DISP</b>     | YES                              | Select SERIAL DISP function                  | NO                 | Reject SERIAL DISP function              |
| 6.7.10 asp. autom. $\Sigma$ VOL ? Y/N | YES                              | Automatic pick up of reagent                 | NO                 | Manual pick up of reagent                |
| 6.7.11 disp automatic ? Y/N           | YES                              | Automatic dispense of reagent                | NO                 | Manual dispense of reagent               |
| 6.7.12 TIME DELAY 1 s                 | YES                              | Memorizes time delay of 1 sec.               | NO                 | Select any time between 0.1 and 9.9 sec. |

| Display                                         | Action       | Result                                 | Alternative Action | Alternative Result Comment      |
|-------------------------------------------------|--------------|----------------------------------------|--------------------|---------------------------------|
| 6.7.13 REAGENT # 1 $\mu\text{L}$                | enter volume | Volume displayed                       | NO                 | Terminates reagent volume entry |
| 6.7.14 REAGENT # 1 ( $V_{R1}$ ) $\mu\text{L}$   | YES          | Volume memorized                       | NO                 | No effect                       |
| 6.7.15 REAGENT # 2 $\mu\text{L}$                | enter volume | Volume displayed                       | NO                 | Terminates reagent volume entry |
| 6.7.16 REAGENT # 2 ( $V_{R2}$ ) $\mu\text{L}$   | YES          | Volume memorized                       | NO                 | VERIFY ?                        |
| 6.7.17 REAGENT # 10 $\mu\text{L}$               | enter volume | Volume displayed                       | NO                 | Terminates reagent volume entry |
| 6.7.18 REAGENT # 10 ( $V_{R10}$ ) $\mu\text{L}$ | YES          | Volume memorized                       | NO                 | No effect                       |
| 6.7.19 VERIFY ?                                 | YES          | Opportunity to verify or amend program | NO                 | Selects RUN channel             |
| 6.7.20 <b>FUNCT</b> VOL SYR SPEED               | YES          | Shows mode of operation                | NO                 | Selects VOL channel             |
| 6.7.21 SERIAL DISPENSER                         | YES          | Memorizes mode of operation            | NO                 | Selects DIL DDL DISP DIST       |
| 6.7.22 asp. automatic $\Sigma$ VOL              | YES          | Memorizes automatic pick up            | NO                 | Selects manual operation        |
| 6.7.23 dispense automatic                       | YES          | Memorizes automatic dispense           | NO                 | Selects manual operation        |

| Display                                   | Action | Result                                      | Alternative Action | Alternative Result Comment               |
|-------------------------------------------|--------|---------------------------------------------|--------------------|------------------------------------------|
| 6.7.24 TIME DELAY 1 s                     | YES    | Memorizes time delay of 1 s                 | NO                 | Select any time between 0.1 and 9.9 sec. |
| 6.7.25 FUNCT <b>VOL</b> SYR SPEED         | YES    | Displays programmed volumes                 | NO                 | Selects <b>SYR</b> channel               |
| 6.7.26 REAGENT # 1 ( $V_{R1}$ ) $\mu$ L   | YES    | Memorizes ( $V_{R1}$ ) volume               | NO                 | Change # 1 reagent volume                |
| 6.7.27 REAGENT # 2 ( $V_{R2}$ ) $\mu$ L   | YES    | Memorizes ( $V_{R2}$ ) volume               | NO                 | Change # 2 reagent volume                |
| .                                         |        |                                             |                    |                                          |
| 6.7.28 REAGENT # 10 ( $V_{R10}$ ) $\mu$ L | YES    | Memorizes ( $V_{R10}$ ) volume              | NO                 | Change # 10 reagent volume               |
| 6.7.29 FUNCT VOL <b>SYR</b> SPEED         | YES    | Shows microprocessor selected syringe sizes | NO                 | Rejects syringe verification             |
| 6.7.30 A-SYRINGE AD ( $V_1$ ) $\mu$ L     | YES    | A-Syringe value accepted                    | NO                 | Shows next suitable syringe              |
| 6.7.31 MICRO(MACRO)-VALVE                 | YES    | Valve type accepted                         | NO                 | Shows alternative valve type             |
| 6.7.32 FUNCT VOL SYR <b>SPEED</b>         | YES    | Shows computed speeds                       | NO                 | VERIFY ?                                 |
| 6.7.33 asp REAGENT $s_3$ sec              | YES    | Accept computed speed to aspirate reagent   | NO                 | Change speed                             |
| 6.7.34 disp REAGENT $s_4$ sec             | YES    | Accept computed speed to dispense reagent   | NO                 | idem                                     |
| 6.7.35 VERIFY ?                           | NO     | All parameters accepted                     | NO                 | Verification required                    |

| Display                                                             | Action                             | Result                                                      | Alternative Action | Alternative Result Comment                                               |
|---------------------------------------------------------------------|------------------------------------|-------------------------------------------------------------|--------------------|--------------------------------------------------------------------------|
| 6.7.36 <b>RUN</b> OR PROG OR AID                                    | YES                                | Select RUN channel to execute method programmed             | NO                 | Selects PROG channel                                                     |
| 6.7.37 Method # (n)                                                 | YES                                | Execute method (n) just programmed                          | NO                 | Select any other method                                                  |
| 6.7.38 (n) SDISP aut # 1 ( $V_{R1}$ ) $\mu$ L                       | YES                                | Displays method and volumes                                 | NO                 | <b>RUN</b> OR PROG OR AID                                                |
| 6.7.39 START TO CHANGE SYR                                          | Press start button                 | Lower drives to allow syringe fitting                       | Comment            | If syringes already installed are o.k., program continues at 6.7.43      |
| 6.7.40 INST MICRO(MACRO)-VALVE                                      | Install valve and indicate by YES  | Microprocessor memorizes which valve fitted                 | Comment/NO         | Only ever change valve when indicated to do so/NO - continue at 6.7.31   |
| 6.7.41 INST ASYR AD ( $V_1$ ) $\mu$ L                               | Install A-Syringe indicate by YES  | Microprocessor memorizes which A-syringe fitted             | Comment/NO         | Only ever change syringe when indicated to do so/NO - continue at 6.7.30 |
| 6.7.42 INST ANY B-SYRINGE                                           | Install any B-Syr. indicate by YES | Microprocessor memorizes which B-syringe fitted             | Comment/NO         | Only ever change syringe when indicated to do so                         |
| 6.7.43 <b>START</b> $\rightarrow$ zero set                          | Press start button                 | MICROLAB 1000 finds zero position of syringes               |                    |                                                                          |
| 6.7.44 <b>START</b> $\rightarrow$ prime <b>N</b> $\rightarrow$ cont | Press start button                 | A-syringe cycles continuously until start button is pressed | NO                 | No prime cycle required                                                  |

| Display                                | Action              | Result                                      | Alternative Action | Alternative Result Comment                |
|----------------------------------------|---------------------|---------------------------------------------|--------------------|-------------------------------------------|
| 6.7.45 endprime⇒PRESS START            | Press start button  | Prime cycle stops                           |                    |                                           |
| 6.7.46 asp ( $V_{RT}$ ) $\mu$ L        | automatically       | A-Syringe aspirates ( $V_{RT}$ ) volume     | YES/NO             | Enters stop routine                       |
| 6.7.47 disp # 1 ( $V_{R1}$ ) $\mu$ L   | automatically       | A-Syringe dispenses ( $V_{R1}$ ) volume     | YES/NO             | idem.                                     |
| 6.7.48 disp # 2 ( $V_{R2}$ ) $\mu$ L   | automatically       | A-Syringe dispenses ( $V_{R2}$ ) volume     | YES/NO             | idem                                      |
| .                                      |                     |                                             |                    |                                           |
| .                                      |                     |                                             |                    |                                           |
| .                                      |                     |                                             |                    |                                           |
| 6.7.49 disp # 10 ( $V_{R10}$ ) $\mu$ L | automatically       | A-Syringe dispenses ( $V_{R10}$ ) volume    | YES/NO             | idem.                                     |
| 6.7.50 As 6.7.46                       | Continues operation |                                             | Comment            | Operation repeats until YES/NO is touched |
| 6.7.51                                 | YES/NO              | Stops routine                               |                    |                                           |
| 6.7.52 <b>STOP</b> OR CONTINUE         | YES                 | SERIAL DISP terminated                      | NO                 | Stops routine                             |
| 6.7.53 <b>RUN</b> OR PROG OR AID       |                     | MICROLAB 1000 ready to function in any mode |                    |                                           |



## 6.8 Program as Multiple pick up

| Display                            | Action                                    | Result                                       | Alternative Action | Alternative Result Comment                           |
|------------------------------------|-------------------------------------------|----------------------------------------------|--------------------|------------------------------------------------------|
| SELF TEST COMPLETED                | Switch on power                           |                                              | Comment            | BRIEFLY DISPLAYED                                    |
| 6.8.1 <b>RUN</b> OR PROG OR AID    | NO                                        | Reject RUN channel                           | YES                | RUN channel selected                                 |
| 6.8.2 RUN OR <b>PROG</b> OR AID    | YES                                       | PROG channel selected                        | NO                 | Reject PROG channel                                  |
| 6.8.3 METHOD #                     | enter unoccupied location number          | Number displayed                             |                    |                                                      |
| 6.8.4 METHOD (n)                   | YES                                       | Location referenced for method being entered | NO                 | No effect                                            |
| 6.8.5 <b>DIL</b> DDL DISP DIST     | NO                                        | Reject DIL function                          | YES                | Select DIL function                                  |
| 6.8.6 DIL <b>DDL</b> DISP DIST     | NO                                        | Reject DDL function                          | YES                | Select DDL function                                  |
| 6.8.7 DIL DDL <b>DISP</b> DIST     | NO                                        | Reject DISP function                         | YES                | Select DISP function                                 |
| 6.8.8 DIL DDL DISP <b>DIST</b>     | YES                                       | Select DIST function                         | NO                 | As 6.8.5 above                                       |
| 6.8.9 <b>SERUM DIST</b> MULTI ASP  | NO                                        | Reject SERUM DIST function                   | YES                | Select SERUM DIST function                           |
| 6.8.10 SERUM DIST <b>MULTI ASP</b> | YES                                       | Select MULTI ASP function                    | NO                 | As 6.8.9 above                                       |
| 6.8.11 ASP # 1 $\mu$ L             | enter first sample volume to be aspirated | Display first sample volume                  | NO                 | Terminates sample volume entry. Selects wash volume. |

| Display                                     | Action                                   | Result                                 | Alternative Action | Alternative Result Comment                           |
|---------------------------------------------|------------------------------------------|----------------------------------------|--------------------|------------------------------------------------------|
| 6.8.12 ASP # 1 ( $V_{S1}$ ) $\mu\text{L}$   | YES                                      | Memorizes first sample volume          | NO                 | No effect                                            |
| 6.8.13 ASP # 2 $\mu\text{L}$                | enter 2nd sample volume to be aspirated. | Display 2nd sample volume              | NO                 | Terminates sample volume entry. Selects wash volume. |
| 6.8.14 ASP # 2 ( $V_{S2}$ ) $\mu\text{L}$   | YES                                      | Memorizes 2nd sample volume            | NO                 | No effect.                                           |
| 6.8.15 ASP # 10 $\mu\text{L}$               | enter 10th sample volume to be aspirated | Display 10th sample volume             | NO                 | Terminates sample volume entry. Selects wash volume. |
| 6.8.16 ASP # 10 ( $V_{S10}$ ) $\mu\text{L}$ | YES                                      | Memorizes 10th sample volume           | NO                 | No effect.                                           |
| 6.8.17 WASH VOLUME $\mu\text{L}$            | enter wash volume                        | Display wash volume                    | NO                 | No effect.                                           |
| 6.8.18 WASH VOLUME ( $V_W$ ) $\mu\text{L}$  | YES                                      | Memorizes wash volume                  | NO                 | No effect.                                           |
| 6.8.19 VERIFY ?                             | YES                                      | Opportunity to verify or amend program | NO                 | Selects RUN channel                                  |
| 6.8.20 <b>FUNCT</b> VOL SYR SPEED           | YES                                      | Shows mode of operation                | NO                 | Selects VOL channel                                  |
| 6.8.21 MULTIPLE ASPIRATION                  | YES                                      | Memorizes function programmed          | NO                 | As 6.8.5 above.                                      |
| 6.8.22 <b>FUNCT</b> <b>VOL</b> SYR SPEED    | YES                                      | Shows mode of operation                | NO                 | Selects SYR channel                                  |
| 6.8.23 ASP # 1 ( $V_{S1}$ ) $\mu\text{L}$   | YES                                      | Memorizes 1st sample volume            | NO                 | Alter 1st sample volume.                             |

| Display                                                                                                                | Action | Result                                      | Alternative Action | Alternative Result Comment        |
|------------------------------------------------------------------------------------------------------------------------|--------|---------------------------------------------|--------------------|-----------------------------------|
| 6.8.23 ASP # 2 ( $V_{S2}$ ) $\mu\text{L}$                                                                              | YES    | Memorizes 2nd sample volume                 | NO                 | Alter 2nd sample volume.          |
| 6.8.24 ASP # 10 ( $V_{S10}$ ) $\mu\text{L}$                                                                            | YES    | Memorizes 10th sample volume                | NO                 | Alter 10th sample volume.         |
| 6.8.25 WASH VOLUME ( $V_W$ ) $\mu\text{L}$                                                                             | YES    | Memorizes wash volume                       | NO                 | Alter wash volume.                |
| 6.8.26 FUNCT VOL <span style="border: 1px solid black;">SYR</span> <span style="border: 1px solid black;">SPEED</span> | YES    | Shows microprocessor selected syringe sizes | NO                 | Rejects syringe verification.     |
| 6.8.27 A-SYRINGE AD ( $V_1$ ) mL                                                                                       | YES    | A-Syringe value accepted                    | NO                 | Shows next suitable syringe size. |
| 6.8.28 B-SYRINGE DAD ( $V_2$ ) $\mu\text{L}$                                                                           | YES    | B-Syringe value accepted                    | NO                 | Shows next suitable syringe size. |
| 6.8.29 MICRO(MACRO)-VALVE                                                                                              | YES    | Valve type accepted                         | NO                 | Shows alternative valve.          |
| 6.8.30 FUNCT VOL SYR <span style="border: 1px solid black;">SPEED</span>                                               | YES    | Shows computed speeds                       | NO                 | VERIFY ?                          |
| 6.8.31 asp SAMPLE ( $s_1$ ) sec                                                                                        | YES    | Accept computed time to aspirate sample.    | NO                 | Change computed speed             |
| 6.8.32 disp SAMPLE ( $s_2$ ) sec                                                                                       | YES    | Accept computed time to dispense sample.    | NO                 | idem.                             |
| 6.8.33 asp WASH ( $s_3$ ) sec                                                                                          | YES    | Accept computed time to aspirate wash       | NO                 | idem.                             |

| Display                                                             | Action                            | Result                                        | Alternative Action | Alternative Result Comment                                               |
|---------------------------------------------------------------------|-----------------------------------|-----------------------------------------------|--------------------|--------------------------------------------------------------------------|
| 6.8.34 disp WASH ( $s_4$ ) sec                                      | YES                               | Accept computed time to dispense wash         | NO                 | idem.                                                                    |
| 6.8.35 VERIFY ?                                                     | NO                                | All parameters accepted                       | YES                | As 6.8.19 above.                                                         |
| 6.8.36 <b>RUN</b> OR PROG OR AID                                    | YES                               | RUN channel selected                          | NO                 | As 6.8.1 above                                                           |
| 6.8.37 METHOD # (n)                                                 | YES                               | Accepts method (n) just programmed            | NO                 | Select any other method.                                                 |
| 6.8.38 (n) MULTI ASP # 1 ( $V_{S1}$ ) $\mu$ L                       | YES                               | Displays method and volumes                   | NO                 | As 6.8.1 above                                                           |
| 6.8.39 START TO CHANGE SYR                                          | Press start button                |                                               | Comment            | If correct syringes are already installed, program continues at 6.8.43.  |
| 6.8.40 INST MICRO(MACRO)-VALVE                                      | Install valve indicate by YES     | Computer memorizes which valve fitted         | Comment/NO         | Only ever change valve when indicated to do so/NO - continue at 6.8.29   |
| 6.8.41 INST ASYR AD                                                 | Install A-syringe indicate by YES | Computer memorizes which A-syringe fitted     | Comment/NO         | Only ever change syringe when indicated to do so/NO - continue at 6.8.27 |
| 6.8.42 INST BSYR DAD                                                | Install B-syringe indicate by YES | Computer memorizes which B-syringe fitted     | Comment/NO         | Only ever change syringe when indicated to do so/NO - continue at 6.8.28 |
| 6.8.43 <b>START</b> $\rightarrow$ zero set                          | Press start button                | MICROLAB 1000 finds zero position of syringes | NO                 | Press start button.                                                      |
| 6.8.44 <b>START</b> $\rightarrow$ prime <b>N</b> $\rightarrow$ cont | NO                                | Priming of system not required                | YES                | Priming of system required.                                              |
| 6.8.45 asp ( $V_W$ ) $\mu$ L                                        | Press start button                | B-Syringe aspirates wash                      | YES/NO             | STOP OR CONTINUE                                                         |

| Display                                     | Action             | Result                                      | Alternative Action | Alternative Result Comment |
|---------------------------------------------|--------------------|---------------------------------------------|--------------------|----------------------------|
| 6.8.46 asp # 1 ( $V_{S1}$ ) $\mu\text{L}$   | Press start button | A-Syringe aspirates 1st sample              | YES/NO             | idem.                      |
| 6.8.47 asp # 2 ( $V_{S2}$ ) $\mu\text{L}$   | Press start button | A-Syringe aspirates 2nd sample              | YES/NO             | idem.                      |
| .                                           |                    |                                             |                    |                            |
| .                                           |                    |                                             |                    |                            |
| .                                           |                    |                                             |                    |                            |
| 6.8.48 asp # 10 ( $V_{S10}$ ) $\mu\text{L}$ | Press start button | A-Syringe aspirates 10th sample             | YES/NO             | idem.                      |
| 6.8.49 disp ( $V_T$ ) $\mu\text{L}$         | Press start button | A-Syringe dispenses total sample volume     | YES/NO             | idem.                      |
| 6.8.50 disp WASH                            | Press start button | B-Syringe dispenses wash volume             | YES/NO             | idem.                      |
| 6.8.51 As 6.8.45 above                      | Press start button | Cycle repeats                               | YES/NO             | idem.                      |
| 6.8.52                                      | YES/NO             | Enters stop routine                         |                    |                            |
| 6.8.53 <b>STOP</b> OR CONTINUE              | YES                | Terminates MULTI ASP function               | NO                 | Program continues.         |
| 6.8.54 <b>RUN</b> OR PROG OR AID            |                    | MICROLAB 1000 ready to function in any mode |                    |                            |

# 6.9 To find a stored method or an unoccupied location

| Display                         | Action | Result                   | Alternative Action | Alternative Result Comment |
|---------------------------------|--------|--------------------------|--------------------|----------------------------|
| 6.9.1 <b>RUN</b> OR PROG OR AID | NO     | Reject RUN channel       | YES                | Select RUN channel         |
| 6.9.2 RUN OR <b>PROG</b> OR AID | NO     | Reject PROG channel      | YES                | Select PROG channel        |
| 6.9.3 RUN OR PROG OR <b>AID</b> | YES    | Select AID channel       | NO                 | Reject AID channel         |
| 6.9.4 <b>REVIEW</b> OR DELETE   | YES    | Select REVIEW feature    | NO                 | Reject REVIEW feature      |
| 6.9.5 (overview of location 1)  | YES    | Advances location number | NO                 | Returns programme to 6.9.1 |
| 6.9.6 (overview of location 2)  | YES    | Advances location number | NO                 | idem                       |
| 6.9.7 (overview of location 50) | YES    | Advances location number | NO                 | idem                       |
| 6.9.8 As 6.9.1                  | YES    | Selects RUN channel      | NO                 | Rejects RUN channel        |

# 6.10 To delete one or all programmes

| Display                          | Action                               | Result                                                   | Alternative Action | Alternative Result Comment                    |
|----------------------------------|--------------------------------------|----------------------------------------------------------|--------------------|-----------------------------------------------|
| 6.10.1 <b>RUN</b> OR PROG OR AID | NO                                   | Reject RUN channel                                       | YES                | Select RUN channel                            |
| 6.10.2 RUN OR <b>PROG</b> OR AID | NO                                   | Reject PROG channel                                      | YES                | Select PROG channel                           |
| 6.10.3 RUN OR PROG OR <b>AID</b> | YES                                  | Select AID channel                                       | NO                 | Reject AID channel                            |
| 6.10.4 <b>REVIEW</b> OR DELETE   | NO                                   | Reject REVIEW feature                                    | YES                | Select REVIEW feature                         |
| 6.10.5 REVIEW OR <b>DELETE</b>   | YES                                  | Select DELETE feature                                    | NO                 | Reject DELETE feature                         |
| 6.10.6 DEL ALL METHOD ? Y/N      | NO                                   |                                                          |                    |                                               |
| 6.10.7 METHOD #                  | Enter number of method to be deleted | Number displayed                                         | YES                | FOR DEL PRESS "."<br>(deletes all programmes) |
| 6.10.8 METHOD (n)                | YES                                  | Computer eliminates entire method from location selected |                    |                                               |
| 6.10.9 <b>RUN</b> OR PROG OR AID | YES                                  | Select RUN channel                                       | NO                 | Reject RUN channel                            |



## SECTION 7

### Fitting the valve block

- 7.1.1 Two types of valve block can be used with the MICROLAB 1000, designated MICRO and MACRO. They should only be fitted or exchanged when indicated by the microprocessor.
- 7.1.2 Lift display cover to give access to the valve block and syringe connections. Insert the two white guides into the two valve drive holes situated beneath the display with the two syringe coupling threaded holes downward. It is not necessary to align the drives with the valves, this will automatically take place during operation. Secure by inserting the central locking screw.

To remove the valve block see section 7.3.1.

### Fitting the syringes

- 7.2.1 Great care must be taken when fitting the syringes not to subject them to any side stress.

Syringes should only be fitted or exchanged when indicated by the microprocessor. No attempt should be made to move the syringe drives by hand, they are programmed to be at their lowest position whenever it is necessary to change or fit different syringes.

- 7.2.2 The MICROLAB 1000 is designed to be used with AccuDi® syringes, a reagent syringe is fitted to the "A" drive (left hand side) and a "D" type, sample syringe to the "B" (right hand side). Sizes to be fitted will be indicated by the microprocessor during dialogue routine. A mechanical stop is fitted to the plunger of the sample syringe to ensure accurate alignment with its side connector outlet in the zero position.
- 7.2.3 Lift display cover to give access to valve block and syringe connections. Loosen the left hand plunger locking ring. Take a reagent/diluent syringe of the volume indicated on the display and pull out the plunger about 2 cm. Holding the syringe vertical with one hand, insert the plunger button into the plunger locking ring assembly and lightly tighten with the other. Extend the syringe and carefully fit barrel into the valve block by twisting with a clockwise motion until tight. Finally tighten the plunger locking ring firmly by hand.

The sample syringe is fitted in a similar manner, but in this case the threaded portion is free to rotate on the barrel, this ensures that the side connector can always be set to face toward the right side of the instrument.

To remove syringes follow the fitting instructions in the reverse order, they may only be dismantled when indicated by the microprocessor.

### To exchange the valve block

- 7.3.1 Dismount syringes, see section 7.2.3. Remove the locking screw from the centre of the valve block and pull the block horizontally away from the main instrument.

To mount another valve block follow 7.1.2.

### **To exchange syringes**

- 7.4.1 The syringes may only be exchanged when indicated by the microprocessor. Dismount syringes fitted and mount the alternatives in accordance with section 7.2.3.

### **To replace defective syringes**

- 7.5.1 If the instrument is running, switch off the mains power and allow 5 seconds to elapse before continuing.

Switch on the power and after the SELF CHECK the MICROLAB 1000 will show the RUN OR PROG OR AID display. Select RUN by touching YES and also the METHOD #(n) used by YES again. Now the display shows START zero set. Press any key but not the start button and MICROLAB 1000 will enter the START TO CHANGE SYR routine. Proceed by following the display dialogue in the normal way, see section 7.2.3. If the syringe size as indicated on the display is not available, touch NO. The instrument will enter the programme mode and will allow a change of syringe size, however, some performance will be lost.

Follow the dialogue in the usual way to program the specific method required.

It may be of help to read the appropriate function under section 6.

| FUNCTION          | SECTION |
|-------------------|---------|
| Diluter           | 6.2.36  |
| Double Diluter    | 6.3.40  |
| Dispenser         | 6.4.31  |
| Serum Distributor | 6.5.41  |
| Serial Diluter    | 6.6.44  |
| Serial Dispenser  | 6.7.39  |
| Multiple pick up  | 6.8.39  |

### **Aspiration tube**

- 7.6.1 The aspiration tubing is fitted directly to the left side of the valve block, using the flared M6 screw connector and tightening firmly by hand. A 2 mm x 1 mm x 500 mm is supplied but this may be shortened if expensive or small volume of reagent is to be used. See also section 7.7.1.

### **Sample/delivery tube**

- 7.6.2 The sample tubing is connected to the side outlet of the sample syringe, using the flared M6 screw connector. The tapered delivery tip should be threaded through handactuator and held in place by lightly tightening the end section. The tube and the cable may be held together using the small clips provided. See also section 7.7.1.

## Selection of tubing

- 7.7.1 If volumes greater than 5 mL are contemplated the aspiration and delivery tubes should be replaced with accessory number 240000/240360 tubing. This has a larger internal diameter than the standard tubing supplied and will produce less resistance to the flow. An alternative would be to reduce the speed of aspiration and delivery by programming. However when it is required to use sample volumes in excess of 700  $\mu$ L the larger delivery tube, number 240360 should be used to avoid the sample entering the syringe, as this may cause contamination producing unacceptable carry-over.

## Cleaning of the FEP tubing

- 7.7.2 After using the MICROLAB 1000 for a longer period, the innerwall contamination of the dispensing tubing installed in the handactuator may disturb the clean separation of the programmed air gap between the sample and reagent. The entire liquid system shall then be cleaned, by flushing several times with a fresh solution of 8 - 13 % sodium hypochlorite, using the PRIME cycle. After cleaning flush system at least 5 times with distilled water before using again with reagents.

### Important:

This procedure should be done at least once a week.

## SECTION 8

### Technical Specifications

|                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                          |          |                          |          |                         |            |                         |          |                         |          |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|----------|--------------------------|----------|-------------------------|------------|-------------------------|----------|-------------------------|----------|
| Methods of operation             | Dispenser, Multiple Dispenser, Diluter, Double Diluter, Serial Diluter, Serum Distributor, Multiple Pick up.                                                                                                                                                                                                                                                                                                                                                                                            |                          |          |                          |          |                         |            |                         |          |                         |          |
| Syringe Drives                   | Each syringe is driven by a high resolution stepper motor and a precision lead screw. An auxiliary synchronmotor drives the valve.                                                                                                                                                                                                                                                                                                                                                                      |                          |          |                          |          |                         |            |                         |          |                         |          |
| Syringe sizes                    | <p>Sample syringe:<br/>50, 100, 250, 500 <math>\mu</math>L, 1.0, 2.5 and 5 mL</p> <p>Diluent syringe:<br/>250, 500 <math>\mu</math>L, 1.0, 2.5, 5.0, 10.0 and 25.0 mL</p> <p>Optimal syringe sizes are calculated by the microprocessor, but the operator may alter this selection.</p>                                                                                                                                                                                                                 |                          |          |                          |          |                         |            |                         |          |                         |          |
| Valves                           | <p>MICRO or MACRO valve</p> <p>Selected by the microprocessor, but may be altered by the user.</p>                                                                                                                                                                                                                                                                                                                                                                                                      |                          |          |                          |          |                         |            |                         |          |                         |          |
| Tubing                           | <p>Standard aspiration tubing:<br/>2 x 1 x 500 mm, (dead volume approx. 390 <math>\mu</math>L)</p> <p>Standard probe tubing:<br/>2 x 1 x 900 mm, (dead volume approx. 700 <math>\mu</math>L)</p> <p>Macro aspiration tubing:<br/>3 x 2 x 500 mm, (dead volume approx. 1.5 mL)</p> <p>Macro probe tubing:<br/>3 x 2 x 900 mm, (dead volume approx. 2.5 mL)</p>                                                                                                                                           |                          |          |                          |          |                         |            |                         |          |                         |          |
| Method storage                   | Maximum 50 methods. These are not lost when power is switched off and may be recalled for use at any time. Memories are protected over 8 years.                                                                                                                                                                                                                                                                                                                                                         |                          |          |                          |          |                         |            |                         |          |                         |          |
| Aspirating and dispensing speeds | Range from 1 to 20 seconds for full stroke of syringe plunger. Following speeds are selectable: 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8, 10, 12, 14, 16, 18 and 20 secs. The requested speeds can be programmed individually for the aspiration of reagent and sample, as well as for the dispensing. Optimum speed is computed by microprocessor, but operator may alter this selection by directly keying in modified values.                                                                                |                          |          |                          |          |                         |            |                         |          |                         |          |
| Display                          | 20 Alpha-numeric characters. Used in dialogue with the operator, i.e. display messages, instructions and positions in program.                                                                                                                                                                                                                                                                                                                                                                          |                          |          |                          |          |                         |            |                         |          |                         |          |
| Accuracy                         | <p>Method of operation: Diluter/Dispenser, applicable to syringes <math>&gt; 250 \mu</math>L</p> <table> <tr> <td>30 % of syringe capacity</td><td><math>&lt; 1 \%</math></td></tr> <tr> <td>10 % of syringe capacity</td><td><math>&lt; 1 \%</math></td></tr> <tr> <td>5 % of syringe capacity</td><td><math>&lt; 1.2 \%</math></td></tr> <tr> <td>2 % of syringe capacity</td><td><math>&lt; 2 \%</math></td></tr> <tr> <td>1 % of syringe capacity</td><td><math>&lt; 3 \%</math></td></tr> </table> | 30 % of syringe capacity | $< 1 \%$ | 10 % of syringe capacity | $< 1 \%$ | 5 % of syringe capacity | $< 1.2 \%$ | 2 % of syringe capacity | $< 2 \%$ | 1 % of syringe capacity | $< 3 \%$ |
| 30 % of syringe capacity         | $< 1 \%$                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                          |          |                          |          |                         |            |                         |          |                         |          |
| 10 % of syringe capacity         | $< 1 \%$                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                          |          |                          |          |                         |            |                         |          |                         |          |
| 5 % of syringe capacity          | $< 1.2 \%$                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                          |          |                          |          |                         |            |                         |          |                         |          |
| 2 % of syringe capacity          | $< 2 \%$                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                          |          |                          |          |                         |            |                         |          |                         |          |
| 1 % of syringe capacity          | $< 3 \%$                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                          |          |                          |          |                         |            |                         |          |                         |          |

|           |                                                                              |
|-----------|------------------------------------------------------------------------------|
| Precision | Method of operation: Diluter/Dispenser, applicable to syringes > 250 $\mu$ L |
|           | 30 % of syringe capacity < 0.2 %                                             |
|           | 10 % of syringe capacity < 0.4 %                                             |
|           | 5 % of syringe capacity < 0.5 %                                              |
|           | 2 % of syringe capacity < 0.8 %                                              |
|           | 1 % of syringe capacity < 1.5 %                                              |

These specifications are statistically based on the test results of more than 500 instruments manufactured and are exceeding the requirements of DIN 12650.

|                   |                                                              |
|-------------------|--------------------------------------------------------------|
| Interface         | RS 232 C (Baud rate 2400)                                    |
| Mains supply      | 100/120/220/240 Volts - 50/60 Hz<br>Voltages are selectable. |
| Power consumption | 40 VA                                                        |
| Fuse ratings      | at 100/120 Volts - 2 Amps<br>at 220/240 Volts - 1 Amp        |
| Dimensions        | 170 x 230 x 300 mm                                           |
| Weight            | 8.5 kg                                                       |

We reserve the right to alter these specifications without notice.

Figure 9.1.1







Figure 9.1.2

Display cover





## SECTION 10

### Description of RS 232 C interface

#### 10.1 Introduction

The MICROLAB 1000 is equipped with a RS 232 C interface which includes transmitted data ( $T_x$ ), received data ( $R_x$ ) and ground (GND).

After switching on the MICROLAB 1000 the display shows RUN OR PROG OR AID. To activate the interface, any ASCII-character has to be sent to the interface before any key on the MICROLAB 1000 keyboard has been touched. The display will then show EXTERNAL. In this mode only single movements can be programmed at a time.

The MICROLAB 1000 may be cascaded with one or more MICROLAB M's. The MICROLAB 1000 has to be the last unit in such a cascade. A MICROLAB unit acts always as a slave.

#### 10.2 Protocol

##### 10.2.1 Line establishment

Transmission is composed by "EOT  $\emptyset$  ENQ" (enquiry) and terminated with "EOT" (end of transmission). Transmission is composed after the MICROLAB 1000 has acknowledged with " $\emptyset$  ACK" (acknowledge). Within the "ENQ" ... "EOT" frame no other unit in a cascade can be addressed.

Example:

MICROLAB 1000 should be activated

|                 |   |                 |
|-----------------|---|-----------------|
| Master          |   | MICROLAB 1000   |
| $\emptyset$ ENQ | → |                 |
|                 | ← | $\emptyset$ ACK |

.

.

.

Command string

.

.

.

EOT →

##### 10.2.2 Communications

Data will be transferred in blocks using 7-bit-ASCII-codes (DIN 66003). The first character of a block is "STX" (start of text). The signal "ETX" (end of text) followed by a BCC (block check character) terminates the block.

Each data block transfer is acknowledged by the receiver:

|                               |                                                  |
|-------------------------------|--------------------------------------------------|
| "ACK" (acknowledge):          | after a successful transfer                      |
| "NAK" (negative acknowledge): | faulty transfer, requests a transfer repetition, |

No answer

in case "ETX" or "BCC" was faulty.

The master must possess a termination criteria, i.e. "Time Out" and a repetition-counter. The MICROLAB 1000 does not possess a termination criteria.

Examples:

successful transfer

|                 |   |               |
|-----------------|---|---------------|
| Master          |   | MICROLAB 1000 |
| STX O G ETX BCC | → |               |
|                 | ← | ACK           |

faulty transfer (transmission error)

|                 |   |               |
|-----------------|---|---------------|
| Master          |   | MICROLAB 1000 |
| STX O G ETX BCC | → |               |
|                 | ← | NAK           |
| STX O G ETX BCC | → | (repetition)  |
|                 | ← | ACK           |

garbled transfer (e.g. ETX faulty)

|                 |   |               |
|-----------------|---|---------------|
| Master          |   | MICROLAB 1000 |
| STX O G ETX BCC | → |               |
|                 |   | no answer     |
| time out        |   |               |
| STX O G ETX BCC | → | (repetition)  |
|                 | ← | ACK           |

### 10.2.3 Data communications-check

To each 7-bit character a parity bit is added. The parity is even. "ETX" (end of text) is followed by a block check character (BCC), which represents the check sum of the block (exclusive of "STX", but inclusive of "ETX"). The parity is odd.

Example:

|     |       |             |            |
|-----|-------|-------------|------------|
|     |       | Even parity |            |
|     |       | ←           |            |
| STX | (02H) | 1 000 0010  |            |
| O   | (4FH) | 1 100 1111  |            |
| G   | (47H) | 0 100 0111  |            |
| ETX | (03H) | 0 000 0011  |            |
| BCC |       | 0 100 0100  | odd parity |

STX has priority and deletes a previous, not terminated string.

#### 10.2.4 Command list

| ASCII | Hex |                                                               |
|-------|-----|---------------------------------------------------------------|
| 'O'   | 4F  | turn valve to <u>O</u> utput                                  |
| 'I'   | 49  | turn valve to <u>I</u> nput                                   |
| 'P'   | 50  | left syringe - <u>P</u> ick up                                |
| 'D'   | 44  | left syringe - <u>D</u> ispense                               |
| 'p'   | 70  | right syringe - <u>p</u> ick up                               |
| 'd'   | 64  | right syringe - <u>d</u> ispense                              |
| 'Z'   | 5A  | initialize                                                    |
| 'S'   | 53  | Status request                                                |
| 'G'   | 47  | <u>G</u> o (Start)                                            |
| 'K'   | 4B  | activate <u>K</u> eyboard / initial status RUN OR PROG OR AID |

# Parameters

| ASCII | HEX |
|-------|-----|
| Ø     | 30  |
| 1     | 31  |
| 2     | 32  |
| 3     | 33  |
| 4     | 34  |
| 5     | 35  |
| 6     | 36  |
| 7     | 37  |
| 8     | 38  |
| 9     | 39  |
| :     | 3A  |
| ;     | 3B  |
| <     | 3C  |
| =     | 3D  |
| >     | 3E  |
| ?     | 3F  |

A syringe command is followed by a parameter block.

## Format:

Ssssv

S syringe command (P,D,p,d)  
s number of steps of stepper motor

$$ssss = \frac{\text{volume } (\mu\text{L})}{\text{syringe capacity } (\mu\text{L})} \times 1000$$

v full stroke = 1000 steps  
speed  
range 30 ... 3F

|             |            |
|-------------|------------|
| Ø = 1 sec   | 8 = 7 sec  |
| 1 = 1.5 sec | 9 = 8 sec  |
| 2 = 2 sec   | : = 10 sec |
| 3 = 2.5 sec | ; = 12 sec |
| 4 = 3 sec   | < = 14 sec |
| 5 = 4 sec   | = = 16 sec |
| 6 = 5 sec   | > = 18 sec |
| 7 = 6 sec   | ? = 20 sec |

indicated time is per full stroke of syringe

## 10.2.5 Status reply

Upon a status request 'S' the MICROLAB 1000 will immediately return a 2 byte status information.

### Note:

No status reply is transmitted in case 'S' was acknowledged by 'NAK'.

# Format:

STX B<sub>1</sub> B<sub>2</sub> ETX BCC

## Byte 1:

Bit

|   |        |                                             |
|---|--------|---------------------------------------------|
| 7 | parity |                                             |
| 6 | 1      |                                             |
| 5 | 'O'    |                                             |
| 4 | 'I'    |                                             |
| 3 | 'P'    | 1 = command running<br>0 = command executed |
| 2 | 'D'    |                                             |
| 1 | 'p'    |                                             |
| 0 | 'd'    |                                             |


## Byte 2:

Bit

|   |        |                          |
|---|--------|--------------------------|
| 7 | parity |                          |
| 6 | 1      |                          |
| 5 | 0      | not used                 |
| 4 | 0      | not used                 |
| 3 | =1     | command 'Z' busy         |
| 2 | =1     | syntax error             |
| 1 | =1     | left syringe overloaded  |
| 0 | =1     | right syringe overloaded |

## Example:

Turn valve to 'INPUT' position



whereby  $\hat{=}$ 

|     |                   |
|-----|-------------------|
| B1  | B2                |
| P @ | 10100000 10000000 |

## 10.2.6 Important notes

- 'Z' has highest priority, it can interrupt any running operation.
- Of the commands O,I,Z,S,K only one can be transmitted at a time.
- P or D can be combined with p or d. Also while one syringe is in process, the second may be released. Except this case, no motion command (O,I,P,D,p,d) is accepted during the busy state of the instrument. A syntax error would be set.
- The execution of the commands O,I,P,D,p,d,Z can be started in 3 different ways:


1. transmitting G together with the command in the same string  
immediate execution
  2. transmitting G in a separate string after the command string  
execution upon receipt of G
  3. no G is transmitted  
execution is released by an external switch (e.g. handactuator)
- After power on and after a syringe overload has occurred no other command than Z will be accepted.
  - The command K returns the control to the keyboard of the MICRO-LAB 1000. This is effective only after termination by EOT.
- Note:  
K acts like a hardware reset and does not pay attention to running commands.

### 10.3 Hardware specifications

The RS 232 C is accessible at the 8 pin DIN socket at the back panel of the MICROLAB 1000.

Pin configuration:

| Pin |                  |              |
|-----|------------------|--------------|
| 1   | LED K            | Handactuator |
| 2   | START SWITCH IN  |              |
| 3   | START SWITCH OUT |              |
| 4   | LED A            | RS 232 C     |
| 5   | GND              |              |
| 6   | Tx               |              |
| 7   | Rx               |              |
| 8   | GND              |              |




An optional Y piece adapter allows the use of the handactuator and the interface at the same time.

Interface specifications:

|              |      |
|--------------|------|
| Baud rate :  | 2400 |
| Parity :     | even |
| Start bit :  | 1    |
| Data bits :  | 7    |
| Parity bit : | 1    |
| Stop bits :  | 2    |

# Branching diagram



- 1) The control station is master station in case no polling has taken place  
R<sub>1</sub>: Recovery after time control is over by means of the control station



## SECTION 11

## Test Certificate

MICROLAB 1000, Serial Number 2988 Date 20.1.88

Results obtained are the average of 15 separate weighings of water at 10 % part volume (22° C)

2500  $\mu$ L A-Syringe, DISP mode: accuracy: -0.11 %  
precision: 0.01 %  
500  $\mu$ L B-Syringe, DIL mode: accuracy: 0.06 %  
precision: 0.26 %

Quality controller: D. Cels-

2500  $\mu$ L A Syringe  
250  $\mu$ L B-Syringe .